

AUG 16, 2023


U54 SCENT 10x Genomic Single-Cell Multiome

Karen Abramson¹, Simon Gregory^{1,2}

¹Duke Molecular Physiology Institute, Durham, NC, USA;

²Department of Neurology, Duke University, Durham, NC, USA

Cellular Senescence Network (SenNet) Method Development Community



valerie.bekker

OPEN ACCESS



ABSTRACT

This document outlines the 10X Genomics Single-Cell Multiome Protocol used for normal lung and colon tissues for SCENT at Duke University.

DOI:

dx.doi.org/10.17504/protocols.io.81wgbxdpolpk/v1

Protocol Citation: Karen Abramson, Simon Gregory 2023. U54 SCENT 10x Genomic Single-Cell Multiome. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.81wgbxdpolpk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jun 07, 2023

10X Genomics Single-Cell Multiome Protocol

- 1** Nuclei were isolated using this 10x Genomics kit - Chromium Nuclei Isolation with RNase Inhibitor Kit, 16rxns, and this [protocol](#)
- 2** Multiome libraries were generated using these kits - Chromium Next GEM Single Cell Multiome ATAC + Gene Expression Reagent Bundle and this [protocol](#)
- 3** Libraries were sequenced on Illumina sequencers to a read depth of at least 50k reads/nucleus for Gene Expression with read lengths of 28x10x10x90 and at least 25k reads/nucleus for ATAC with read lengths of 50x8x24x49.