




JAN 22, 2024

## 1M Potassium Phosphate Buffer, pH 6.0

 In 1 collection

Adrien Assie<sup>1</sup>, Buck Samuel<sup>1</sup>

<sup>1</sup>Baylor College of Medicine

Samuel Lab



Adrien Assie  
Baylor College of Medicine

OPEN  ACCESS



### ABSTRACT

Potassium Phosphate buffer recipe

**Protocol Citation:** Adrien Assie, Buck Samuel 2024. 1M Potassium Phosphate Buffer, pH 6.0. **protocols.io**  
<https://protocols.io/view/1m-potassium-phosphate-buffer-ph-6-0-zyzf7x6>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** Apr 10, 2019

**Last Modified:** Jan 22, 2024

**PROTOCOL integer ID:** 22265

1 Start with  700 mL water

2 108.3 g KH<sub>2</sub>PO<sub>4</sub>

3 35.6 g K<sub>2</sub>HPO<sub>4</sub>

4 Adjust to  1000 mL water

5 Autoclave