

NOV 23, 2023

## Open field

Núria

Miquel Vila<sup>1</sup>, Peñuelas<sup>1</sup>, Ariadna Laguna<sup>1</sup>

<sup>1</sup>Vall d'Hebron Research Institute

Vilalab Public



Miquel Vila

**ABSTRACT** 

Open field





DOI:

dx.doi.org/10.17504/protocol s.io.6qpvr32r2vmk/v1

**Protocol Citation: Miquel** Vila, Núria Peñuelas, Ariadna Laguna 2023. Open field. protocols.io

https://dx.doi.org/10.17504/p rotocols.io.6qpvr32r2vmk/v1

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License. which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status: Working** We use this protocol and it's working

Created: Nov 22, 2023

Last Modified: Nov 23,

2023

Oct 23 2023

## **PROTOCOL integer ID:** 91310

1	Place the mice in a square open field arena of 80 x 80 cm of white methacrylate under brightly lit (300
	lux).

- 2 Track the center and the periphery using a tracking software (SMART 3.0 Panlab-Harvard Apparatus).
- 3 Clean the arena every time between animals.