

DEC 01, 2023

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocol s.io.3byl4q1kzvo5/v1

Protocol Citation: Pranay Srivastava 2023. Chemicals and treatment paradigms. **protocols.io**

https://dx.doi.org/10.17504/protocols.io.3byl4q1kzvo5/v1

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: Nov 30, 2023

Last Modified: Dec 01,

2023

Chemicals and treatment paradigms

Pranay Srivastava¹

¹Massachusetts General Hospital, ASAP

ASAP Collaborative Research Network

Team Chen

1 more workspace ↓



Fang Zhou

Massachusetts General Hospital, Aligning Science Across Park...

ABSTRACT

This protocol is to enhance microglia response and nigral dopaminergic cell death in an acute MPTP model following peripheral inflammation resulting from single i.p. injection of LPS at 2 mg/kg (García-Domínguez et al. 2018).

PROTOCOL integer ID:

91665

Keywords: ASAPCRN

- Mice were randomly divided into MPTP+LPS+NDP-MSH, MPTP+LPS, and control groups to receive i.p. once daily MPTP·HCl (Millipore Sigma, Cat# M0896; 20 mg/kg) or saline and LPS (Millipore Sigma, Cat# L4391; 1 mg/kg) or PBS from day 1 to day 4.
- 2 NDP-MSH (Genscript, Cat# RP10658; 400 μg/kg) or PBS was injected from day 1 to day 12.
- 3 Mice were tested for behavioral activities and were sacrificed thereafter on day 12.

To study the role of Tregs

- **4** Animals were treated with anti-mouse CD25 monoclonal antibody (clone PC61, Biolegend, Cat# 102059; 400 μg/mouse) or isotype control (Biolegend, Cat# 401916; 100μg/mice) for 3 alternate days, 1 week before the start of experiment.
- 5 Mice were subsequently treated with MPTP, LPS and NDP-MSH as described above.
- **6** Another dose of anti-mouse CD25 monoclonal antibody or isotype control was administered 2 days before the sacrifice.