



Sep 03, 2020

# Qbiotix protocol

vito<sup>1</sup><sup>1</sup>Qbiotix

In Development

[dx.doi.org/10.17504/protocols.io.bks3kwgn](https://dx.doi.org/10.17504/protocols.io.bks3kwgn)

vito vito

DOI

[dx.doi.org/10.17504/protocols.io.bks3kwgn](https://dx.doi.org/10.17504/protocols.io.bks3kwgn)

PROTOCOL CITATION

vito 2020. Qbiotix protocol. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bks3kwgn>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 03, 2020

LAST MODIFIED

Sep 03, 2020

PROTOCOL INTEGER ID

41531

- 1- collect saliva vial 1 mL
- 2- open vial and scan barcode
- 3- dispense 250 µl sample to single well in plate with a pipette
- 4- dispense 250 µl of magnetic bead solution in well with a pipette
- 5- incubate at controlled temperature on the plate reader ⌚ 00:10:00
- 6- place plate on magnetic separator ⌚ 00:00:30

7 7- insert plate on analyzer (reading speed depends on chosen plate reader instrument)

8 6- read results