

JUN 01, 2023

OPEN ACCESS

DOI:

dx.doi.org/10.17504/protocol s.io.6qpvre38plmk/v1

Protocol Citation: Amy Gladfelter 2023. CHROMagar Malassezia Medium (CHROM). **protocols.io**

https://dx.doi.org/10.17504/protocols.io.6qpvre38plmk/v1

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: Nov 03, 2018

Last Modified: Jun 01, 2023

PROTOCOL integer ID:

17458

1 CHROMagar Malassezia basal medium 🗸 56.3 g

п гс 2 с

CHROMagar Malassezia Medium (CHROM)

In 1 collection

Amy Gladfelter¹

¹Duke University



Amy Gladfelter

ABSTRACT

Media recipe for growth of Malassezia species.

ATTACHMENTS

Malassezia media_1.pdf

protocols.io | https://dx.doi.org/10.17504/protocols.io.6qpvre38plmk/v1

Note

CHROMagar, Paris, France

- 2 Tween 40 🗸 10 mL
- **3** dH₂O
 ☐ 1 L
- 4 Sterilize by autoclaving.