

Jul 01, 2024

## Optimized Protocol for RNA Extraction from Insect Samples Using TRIzol Reagent

DOI

**[dx.doi.org/10.17504/protocols.io.yxvmepx9g3p/v1](https://dx.doi.org/10.17504/protocols.io.yxvmepx9g3p/v1)**

Monique da Silva Bonjour<sup>1</sup>, Ian de Paula Alves Pinto<sup>1</sup>, Angelo José Rinaldi<sup>1</sup>

<sup>1</sup>Department of Biochemistry and Molecular Biology, Universidade Federal de Viçosa - UFV, BIOAGRO/INCT-IPP, Viçosa-MG, Brazil



Monique da Silva Bonjour

Unversidade Federal de Viçosa

---

OPEN  ACCESS



DOI: **[dx.doi.org/10.17504/protocols.io.yxvmepx9g3p/v1](https://dx.doi.org/10.17504/protocols.io.yxvmepx9g3p/v1)**

**Protocol Citation:** Monique da Silva Bonjour, Ian de Paula Alves Pinto, Angelo José Rinaldi 2024. Optimized Protocol for RNA Extraction from Insect Samples Using TRIzol Reagent. **protocols.io** **<https://dx.doi.org/10.17504/protocols.io.yxvmepx9g3p/v1>**

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** June 28, 2024

**Last Modified:** July 01, 2024

**Protocol Integer ID:** 102584



## Abstract

RNA extraction from insect samples is a crucial step in molecular biology studies aimed at understanding gene expression, functional genomics, and various other applications. This study presents an optimized protocol for RNA extraction from insect tissues using TRIzol reagent. The protocol involves homogenizing insect samples in TRIzol, followed by phase separation with chloroform, RNA precipitation with isopropanol, and washing with ethanol. The final RNA pellet is dissolved in RNase-free water. This method ensures high yield and purity of RNA, suitable for downstream applications such as RT-PCR and RNA sequencing. The efficacy of this protocol was validated through spectrophotometric analysis and agarose gel electrophoresis, demonstrating its reliability and efficiency for extracting RNA from a wide variety of insect species

## Materials

- 1- Insects
- 2- 75% ethanol
- 3- TRIzol (Mixture of guanidine thiocyanate and phenol in a single-phase solution)
- 4- crucible and pistil
- 5- chloroform
- 6- Isopropanol
- 7- Microcentrifuge tubes
- 8- Centrifuge

## Steps

- 1 Cover the crucible with liquid nitrogen to keep the sample cold
- 2 Place the sample (about 50mg) and macerate, keeping it cold. Remove the wings and macerate insect samples as they have a strong exoskeleton and the presence of chitin
- 3 Mix the sample with Trizol in the crucible and transfer to a 1.5 mL microtube
- 4 Add 200  $\mu$ L of ice-cold chloroform for phase separation. Shake the tube vigorously by hand for a full 15 seconds. Let sit at room temperature for 2-3 minutes  
  
**NOTE:** You should see the two mixtures separating almost immediately: pink on the bottom and clear on top
- 5 Add 1mL of Trizol and wait for it to thaw
- 6 Centrifuge at 12000 rpm for 15 min at 4° C  
  
**NOTE:** Any denatured proteins might appear as a white interface. Avoiding this interface, carefully move the clear top liquid to a new clean, labeled tube
- 7 Remove 500  $\mu$ L of supernatant and add 400  $\mu$ L of ice-cold isopropanol to precipitate the RNA (always 100  $\mu$ L less)
- 8 Leave 10 minutes at room temperature  
  
**NOTE:** Leaving to precipitate in the -20°C is not recommended because it might cause contaminants to co-precipitate
- 9 Centrifuge at 12000 rpm for 10 min 4° C  
  
**NOTE:** You may see the pellet as a clear, gel-like ball at the base of the tube, or a small white smear
- 10 Remove the isopropanol and wash with 1 mL of 75% ethanol



- 11 Centrifuge at 9500 rpm for 5 min at 4° C
- 12 Discard the alcohol and wait for the pellet to dry for 5 to 10 minutes at room temperature  
**NOTE:** Remove as much of the ethanol as possible
- 13 Resuspend in 30 µL of RNase-free water depending on the size of the pellet formed

## Protocol references

Avramov, M., Gallo, V., Gross, A. et al. A cost-effective RNA extraction and RT-qPCR approach to detect California serogroup viruses from pooled mosquito samples. *Sci Rep* 14, 2339 (2024). <https://doi.org/10.1038/s41598-024-52534-1>

Han, P., Go, M.K., Chow, J.Y. et al. A high-throughput pipeline for scalable kit-free RNA extraction. *Sci Rep* 11, 23260 (2021). <https://doi.org/10.1038/s41598-021-02742-w>

RNA extraction from Trizol. Protocol for insects.

[https://homepages.eawag.ch/~vorburch/Files/RNA\\_extraction\\_TRIZOL\\_Insects\\_ABD\\_April2016.pdf](https://homepages.eawag.ch/~vorburch/Files/RNA_extraction_TRIZOL_Insects_ABD_April2016.pdf)