





Jul 14, 2022

ID16A measurements on pollen samples

Alexandra Pacureanu¹, Ruxandra Cojocaru^{1,2}, Oonagh Mannix^{1,3}

¹ESRF - The European synchrotron, Grenoble, France;

²Okinawa Institute of Science and Technology Graduate University, Onna, Okinawa 904-0495, Japan;

³Helmholtz-Zentrum Berlin für Materialien und Energie, Hahn-Meitner-Platz 1, 14109 Berlin, Germany

1 Works for me Share

dx.doi.org/10.17504/protocols.io.bpfwmjpe

Ruxandra Cojocaru

ESRF - The European Synchrotron, Okinawa Institute of Scienc...

ABSTRACT

X-ray Phase-contrast nano-tomography and X-ray fluorescence microscopy was performed on modern Pinaceae (pine) pollen at the ID16A nano-imaging beamline of the European Synchrotron (ESRF), in France.

DOI

dx.doi.org/10.17504/protocols.io.bpfwmjpe

PROTOCOL CITATION

Alexandra Pacureanu, Ruxandra Cojocaru, Oonagh Mannix 2022. ID16A measurements on pollen samples . **protocols.io** https://dx.doi.org/10.17504/protocols.io.bpfwmjpe

KEYWORDS

phase-contrast, nanotomography, X-rays, nanoscale, pollen, pine pollen, bisaccate pollen, palynology, materials science

LICENSE

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 06, 2020

LAST MODIFIED

Jul 14, 2022



1

PROTOCOL INTEGER ID

44246

MATERIALS TEXT

Modern Pinaceae (pine) pollen was collected from the region of Grenoble (France) 24 February 2018, registered at the Natural History Museum London here. Pollen grains were cleanly extracted from male cones under a microscope. No chemical pre-processing was done before this experiment.

Experimental setup

- 1 Beamline: European Synchrotron Radiation Facility (ESRF) ID16A
 - 185 m long nano-imaging beamline that provides nano-focused beams.
 - Beam size between 30x30 nm to 400x400 um.
 - Can perform X-ray phase-contrast nano-tomography and X-ray fluorescence microscopy (XRF)
 - Hard X-ray photon energies available: 17.05 keV or 33.6 keV.
 - Cryogenic sample preservation available.

Reference:

http://www.esrf.eu/cms/live/live/en/sites/www/home/UsersAndScience/Experiments/XNP/I D16A/over.html (last accessed on January 18th, 2020)

Nano-imaging Beamline

The beamline is part of the ESRF Upgrade project <u>UPBL4 NINA</u> ☑. The 185 m long beamline provides nano-focused beams for analytical imaging. The nano-imaging beamline addresses problems in biology, biomedicine and nano-technology using X-ray fluorescence microscopy and nano-tomography. It is optimised for ultimate hard X-ray focusing of a beam with a large energy bandwidth at specific energies (17keV or 33.6 keV). Cryogenic sample preservation is available for 2D and 3D imaging.

Disciplines

- Life Sciences
- Medicine
- Materials and Engineering
- Environmental Sciences
- Earth and Planetary Sciences
- Physics

Applications

- Biomedicine
- Environmental sciences
- Energy
- Nanotechnology
- Advanced materials

Energy range 17.0 - 33.6 keV

Beam size

Minimum (H x V)

30 x 30 nm²

Maximum (H x V)
400 x 400 μm²

Techniques

- coherent imaging
- phase-contrast imaging
- near-field and far-field ptychography
- X-ray fluorescence microscopy
- scanning transmission X-ray microscopy
- nanotomography

Screenshot from the beamline ID16A ESRF webpage with beamline specs



2

2 Detectors used:

- Imaging: high-resolution imaging detector lens-coupled to a FReLoN F_E230-84 (4096x4096 pixels, 1.5 um pixel size).
- Fluorescence: two six elements silicon drift diode detectors.

Reference:

http://www.esrf.eu/cms/live/live/en/sites/www/home/UsersAndScience/Experiments/XNP/I D16A/over.html (last accessed on January 18th, 2020)

Detectors

- Fluorescence: two six elements silicon drift diode detectors
- Imaging 1: high-resolution imaging detector lens-coupled to a FReLoN F_K4320 (2048x2048 pixels, 1.1 um pixel size)
- Imaging 2: high-resolution imaging detector lens-coupled to a FReLoN
 F_E230-84 (4096x4096 pixels, 1.5 um pixel size) and to a SVCam-HR16070 (4864x3232 pixels, 1.5 um pixel size)
- · Far-field Ptychography: hybrid detector with four modules MAXIPIX detector

Screenshot from the beamline ID16A ESRF webpage with detector specs

- 3 See diagrams and descriptions of experimental setups for:
 - X-ray phase-contrast nano-tomography at ID16A in: M. Hubert, A. Pacureanu, C. Guilloud, Y. Yang, J. C. da Silva, J. Laurencin, F. Lefebvre-Joud, P. Cloetens, Efficient correction of wavefront inhomogeneities in x-ray holographicnanotomography by random sample displacement. Appl Phys Lett 112, 203704 (2018). URL: https://pubs.acs.org/doi/abs/10.1021/acs.analchem.9b04096
 - X-ray fluorescence microscopy at ID16A in: C. Gramaccioni, Y. Yang, A. Pacureanu, N. Vigano, A. Procopio, P. Valenti, L. Rosa, F. Berlutti, S. Bohic, P. Cloetens, Cryonanoimaging of single human macrophage cells: 3D structural and chemical quantification. Anal Chem 92, 4814-4819 (2020). URL: https://aip.scitation.org/doi/full/10.1063/1.5026462

Mounting the samples

4 Sample details can be found in the materials section. Pollen specimens were glued to a metallic Tungsten tip, cut from a 125 micrometer diameter wire, using commercially available gel Superglue. Eyelashes washed in ethanol and glued to toothpicks were used to manipulate the individual pollen grains. Multiple grains were glued to the top of each tip.

- 5 Experimental conditions:
 - Cryogenic conditions: -160 degrees Celsius.
 - High vacuum: 10⁽⁻⁸⁾ mbar.

Cryogenic conditions were chosen to protect the sample from radiation damage.

Experimental parameters

- 6 Nano-tomography:
 - Hard X-ray photon energy used: 17.05 keV.
 - Voxel size: between 20 and 40 nm, depending on the scan and sample.

7 XRF:

- Excitation energy: 17.05 keV.
- Pixel size: 150 nm2.
- Dwell time: 50 ms.
- Spectra recorded with a silicon drift energy dispersive detector.
- Resulting data analysed with <u>PyMCA</u>.
- Raw data available from the ESRF ICAT here: https://doi.esrf.fr/10.15151/ESRF-DC-339743082

Nano-tomography reconstructions

- Data reconstructed using Octave custom code to obtain phase maps.
 - Tomographic volumes computed using <u>PyHST</u>.
 - Specific parameters can be found in the accompanying files for each reconstructed volume
 - Three-dimensional renderings obtained using ParaView software.

Post-processing

The foam cavity sizes in different regions of the reconstructed samples were determined using an adapted version of the post-processing algorithm form <u>Schladitz et al. 2008</u>, originally used for open aluminium foams.

Reference: K. Schladitz, C. Redenbach, T. Sych, M. Godehardt, Microstructural characterisation of open foams using 3d images. Berichte des Fraunhofer ITWM 148 (2008).

The steps implemented in the <u>FIJI/ImageJ</u> software were the following:

- Select a volume of interest (substack) in the desired region (cappa or sacci).
- Straighten if necessary.
- Reslice in top-down direction (radially to the outer surface).
- Binarisation of the substack.

protocols.io

- Compute a distance transform watershed (gray color map).
- Binarisation of the result using an adjusted threshold.
- Inversion (needed or not, depending on ImageJ implementation).
- Watershed again.
- Use Particle Analysis module to automatically identify and measure the cavities.
- Set scale size (1 pixel to voxel size).
- Measure Particles: only cavities above 50 nm2 (basically excludes single pixel "cavities"), using the options "exclude on edges", "show summary", "display".
- Compute histogram, average, standard deviation of resulting list of cavity sizes.
- 11 Density values can be retrieved from tomographic reconstructions, which provide:

$$\omega = -rac{2\pi}{\lambda}\delta$$

where $\lambda=0.72727A^\circ$ is the photon wavelength for an energy of 17.05 keV.

Using Guinier approximation (Z/Approx 1/2):

$$\delta=1.310^{-6}
ho[g/cm^3]\lambda^2[A^\circ]$$

and

$$\omega = -2\pi*1.310^4
ho[g/cm^3]\lambda[A^\circ]$$

Thus, the density is:

$$ho[g/cm^3] = -rac{10^{-2}\omega[cm^{-1}]}{1.3\lambda[A^\circ]2\pi}$$

This value must be corrected by an offset. For this, we consider a region of known density in the reconstruction (air). Here, vacuoles within the pollen grain had the lowest density. By considering these to be filled with with air, the offset can be determined.