

Nov 24, 2025

LC-MS/MS Untargeted Metabolomics Data Processing

DOI

dx.doi.org/10.17504/protocols.io.bf4wjqxe

Kevin Contrepais¹

¹Stanford University

Human BioMolecular Atlas Program (HuBMAP) Method Development Community

Tech. support email: Jeff.spraggins@vanderbilt.edu



Kevin Contrepais

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bf4wjqxe>

Protocol Citation: Kevin Contrepais 2025. LC-MS/MS Untargeted Metabolomics Data Processing. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bf4wjqxe>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: May 06, 2020

Last Modified: May 13, 2020

Protocol Integer ID: 36726

Keywords: HuBMAP, Metabolomics, Stanford University, ms untargeted metabolomics data processing scope, untargeted metabolomics data, metabolomics data, metabolomics dataset, metabolite annotation, m, spectral database,

Abstract

Scope:

To describe the procedure to process untargeted metabolomics data.

Expected outcome/data:

Metabolomics dataset is processed using commercial software and custom R scripts. Metabolite annotation was performed using an in-house library and spectral databases.

Troubleshooting

- 1 Data from each mode were independently analyzed using Progenesis QI software (v2.3) (Nonlinear Dynamics, Durham, NC).
- 2 Metabolic features from blanks and those that didn't show sufficient linearity upon dilution in QC samples ($r < 0.6$) were discarded.
- 3 Only metabolic features present in $>2/3$ of the samples were kept for further analysis.
- 4 Median normalization was applied to correct for differential starting material quantity.
- 5 Missing values were imputed by drawing from a random distribution of low values in the corresponding sample.
- 6 Data from each mode were merged and metabolites were formally identified by matching fragmentation spectra and retention time to analytical-grade standards when possible or matching experimental MS/MS to fragmentation spectra in publicly available databases.
- 7 Metabolite abundances were reported as spectral counts.
- 8 Peak annotation was first performed by matching experimental m/z , retention time and MS/MS spectra to an in-house library of analytical-grade standards. Remaining peaks were identified by matching experimental m/z and fragmentation spectra to publicly available databases including HMDB (<http://www.hmdb.ca/>), MoNA (<http://mona.fiehnlab.ucdavis.edu/>) and MassBank (<http://www.massbank.jp/>) using the R package 'MetID' (v0.2.0) (PMID: 30944337). Briefly, metabolic feature tables from Progenesis QI were matched to fragmentation spectra with a m/z and a retention time window of ± 15 ppm and ± 30 s (HILIC) and ± 20 s (RPLC), respectively. When multiple MS/MS spectra match a single metabolic feature, all matched MS/MS spectra were used for the identification. Next, MS1 and MS2 pairs were searched against public databases and a similarity score was calculated using the forward dot-product algorithm which takes into account both fragments and intensities (PMID: 24222034). Metabolites were reported if the similarity score was above 0.4. Spectra from metabolites of interest were further manually investigated to confirm identification.