

JUN 12, 2023


U54 SCENT 10X Genomics Visium Protocol

Karen Abramson¹, Simon Gregory^{1,2}

¹Duke Molecular Physiology Institute, Durham, NC, USA;

²Department of Neurology, Duke University, Durham, NC, USA

Cellular Senescence Network (SenNet) Method Development Community



valerie.bekker

OPEN ACCESS

Protocol Citation: Karen Abramson, Simon Gregory 2023. U54 SCENT 10X Genomics Visium Protocol . **protocols.io** <https://protocols.io/view/u54-scent-10x-genomics-visium-protocol-cvbaw2ie>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jun 07, 2023

Last Modified: Jun 12, 2023

PROTOCOL integer ID:
83010

ABSTRACT

This document outlines the 10X Genomics Visium Protocol used for normal lung and colon tissues for SCENT at Duke University.

10X Genomics Visium Protocol

- 1 Once slides received, the following 10x Genomics protocol is followed:

https://cdn.10xgenomics.com/image/upload/v1660261285/support-documents/CG000409_Demonstrated_Protocol_VisiumSpatialFFPE_Deparaffin_H_E_RevC.pdf

https://cdn.10xgenomics.com/image/upload/v1660261285/support-documents/CG000407_VisiumSpatialGeneExpressionforFFPE_UserGuide_RevD.pdf

- 2** 10x Genomics kit is used – Visium (v1) Spatial for FFPE Gene Expression Kit, Human Transcriptome.
- 3** Libraries sequenced on Illumina sequencers to a read depth of a least 25k reads/spot, with read lengths of 28x10x10x50.