



FEB 12, 2024

🌐 Elevated body swing test in rats

eduard.bentea¹, María Sanchiz Calvo¹, Veerle Baekelandt¹

¹Katholieke Universiteit Leuven



eduard.bentea

ABSTRACT

Protocol for performing the elevated body swing test in rats. The elevated body swing evaluates motor asymmetry in body swings when the rat is briefly hanged by the tail, and can be used in animal models with unilateral lesions of the nigrostriatal dopaminergic pathway.

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.bp2l6x58zlqe/v1

Protocol Citation: eduard.bentea, María Sanchiz Calvo, Veerle Baekelandt 2024. Elevated body swing test in rats. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bp2l6x58zlqe/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Feb 05, 2024

Last Modified: Feb 12, 2024



PROTOCOL integer ID: 94696

Keywords: ASAPCRN, motor function, behavioral test, rats

Funders Acknowledgement:
ASAP (Aligning Science Across Parkinson's)

Test

1h 0m 5s

- 1 Bring cages to the behavioral room for at least  01:00:00 prior to the test for habituation 1h
- 2 Place rat individually in a new cage
- 3 Wait until attains neutral position with all four limbs touching the bottom of the cage
- 4 Lift vertically by the base of the tail
- 5 Record the first direction of body deviation away from the vertical axis of 30 degree during a maximum  00:00:05 interval 5s
- 6 Repeat the test five times, with a brief rest in between

