




Sep 01, 2022

# Isolation of natural killer (NK) cells from human blood products

Philippa R Kennedy<sup>1</sup><sup>1</sup>University of Minnesota

1 Works for me

 Share[dx.doi.org/10.17504/protocols.io.81wgbp94yvpk/v1](https://dx.doi.org/10.17504/protocols.io.81wgbp94yvpk/v1) Philippa R Kennedy  
University of Minnesota

## ABSTRACT

Standard isolation procedure for peripheral blood mononuclear cells from human blood leukopheresis products

DOI

[dx.doi.org/10.17504/protocols.io.81wgbp94yvpk/v1](https://dx.doi.org/10.17504/protocols.io.81wgbp94yvpk/v1)

## PROTOCOL CITATION

Philippa R Kennedy 2022. Isolation of natural killer (NK) cells from human blood products. **protocols.io**  
<https://protocols.io/view/isolation-of-natural-killer-nk-cells-from-human-bl-bfcvjiw6>



## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Apr 21, 2020

## LAST MODIFIED

Sep 01, 2022

## PROTOCOL INTEGER ID

35957

- 1 Healthy donor blood products (Leukopaks) were obtained from Memorial Blood Bank (Minneapolis, MN). All samples were de-identified and their use was approved by the University of Minnesota and NMDP institutional review board in accordance with the

Declaration of Helsinki.

- 2 Peripheral blood mononuclear cells (PBMCs) were separated through density gradient filtration (Ficoll-Paque; GE Healthcare).
- 3 For experiments with enriched NK cells, PBMCs were processed fresh using negative magnetic bead separation (EasySep Human NK Cell Enrichment Kit; Cat. No:19055, STEMCELL Technologies).
- 4 Experiments with PBMCs or enriched NK cells were performed in RPMI 1640 supplemented with 10% Fetal Bovine Serum and 1% penicillin/streptomycin. Where specified, culture was supplemented with recombinant human IL-15 (National Cancer Institute).