



Nov 18, 2020

Preparo de meio de LB e LB ágar em placa 11/11/20

Mariana Jacques¹¹Universidade Federal do Rio de Janeiro*In Development*

This protocol is published without a DOI.

Phosphate**Mariana Jacques**

PROTOCOL CITATION

Mariana Jacques 2020. Preparo de meio de LB e LB ágar em placa 11/11/20. **protocols.io**
<https://protocols.io/view/preparo-de-meio-de-lb-e-lb-gar-em-placa-11-11-20-bpuhmt6>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 18, 2020

LAST MODIFIED

Nov 18, 2020

PROTOCOL INTEGER ID

44649

Preparo do LB

1Dissolver 25 g de LB em 1 L de H₂O filtrada

Preparo do LB ágar

2Dissolver 12 g de LB ágar 300 mL de H₂O filtrada**3**

Depois de pronto autoclavar os meios LB e LB em placa, em temperatura- controlada.

4

Ao retirar e o meio atingir uma temperatura ideal () adicionar , Ampicilina.

Colocar os meios em placa e guardar.

