



Upload image

Jul 31, 2020

# 16S Bacteria 338F-516P-805R BSA

Forked from 16S Bacteria 338F-516P-805R

Roey Angel<sup>1</sup>, Eva Petrova<sup>1</sup>, Ana Lara<sup>1</sup><sup>1</sup>Soil and Water Research Infrastructure**1** Works for me [dx.doi.org/10.17504/protocols.io.m7qc9mw](https://dx.doi.org/10.17504/protocols.io.m7qc9mw)

SoWa RI Anaerobic and Molecular Microbiology (public)

Tech. support email: [eva.petrova@bc.cas.cz](mailto:eva.petrova@bc.cas.cz) Eva Petrova  
Soil and Water Research Infrastructure

## ABSTRACT

Universal 16S rRNA Bacteria 338F-516P-805R

BAC338F ACT CCT ACG GGA GGC AG , target *E.coli*: 338-354, Yu et al. (2005), B&BBAC516P\* TGC CAG CAG CCG CGG TAA TA, target *E.coli*: 516-536, 'BAC805R GAC TAC CAG GGT ATC TAA TC , target *E.coli*: 785-805, '

\* Probe must be dual-labelled either with 5'-6-FAM, 3'-BHQ1 or any other valid combination.

## THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Yu, Y, Lee, C, Kim, J & Hwang, S (2005) Group-specific primer and probe sets to detect meth-anogenic communities using quantitative real-time polymerase chain reaction. *Biotechnology and Bioengineering* 89:670–679. url:<http://dx.doi.org/10.1002/bit.20347>

## ATTACHMENTS

[Introduction\\_QPCR\\_Strata](#) [AB\\_rt-QPCRguide.pdf](#)  
[gene.pdf](#)

## DOI

[dx.doi.org/10.17504/protocols.io.m7qc9mw](https://dx.doi.org/10.17504/protocols.io.m7qc9mw)

## PROTOCOL CITATION

Roey Angel, Eva Petrova, Ana Lara 2020. 16S Bacteria 338F-516P-805R BSA. **protocols.io**  
[dx.doi.org/10.17504/protocols.io.m7qc9mw](https://dx.doi.org/10.17504/protocols.io.m7qc9mw)

## MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Yu, Y, Lee, C, Kim, J & Hwang, S (2005) Group-specific primer and probe sets to detect meth-anogenic communities using quantitative real-time polymerase chain reaction. *Biotechnology and Bioengineering* 89:670–679. url:<http://dx.doi.org/10.1002/bit.20347>

## FORK FROM

Forked from 16S Bacteria 338F-516P-805R, Eva Petrova

## LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Feb 12, 2018

LAST MODIFIED

Jul 31, 2020

PROTOCOL INTEGER ID

10192

PARENT PROTOCOLS

In steps of

[RNA-Stable Isotope Probing](#)

MATERIALS

NAME	CATALOG #	VENDOR
iQ™ SYBR® Green Supermix	1708880	BioRad Sciences

#### QPCR mixture

1	Reagent	Final conc.	1 tube (20µl)	plate (20µl x 100)
	PCR H <sub>2</sub> O		4.6	460
	iQ™ Supermix	1x	10	1000
	MgCl <sub>2</sub> (25mM)	4.0 mM	0.8*	80*
	BSA (20µg/µl)	0.2 µg/µl	0.2	20
	<b>338F</b> (10µM)	0.5 µM	1.0	100
	<b>805R</b> (10µM)	0.5 µM	1.0	100
	<b>516P</b> (10µM) <sup>†</sup>	0.2 µM	0.4	40
	Template		2	2 x 100

\* Buffer contains MgCl<sub>2</sub> at final conc. of 3.0 mM

#### Program

- 2
1. 95°C – 5'
2. **x 40** {
  - a. 95°C – 30"
  - b. 22°C – 60" take a snapshot}