



16S_PCR_Sequence_Analysis

COMMENTS 0

DOI

dx.doi.org/10.17504/protocols.io.6qpvr49e3gmk/v1

Kenneth Acosta¹

¹Rutgers University

duckweed-bacteria interactions

WORKS FOR ME 1



Kenneth Acosta

ABSTRACT

Protocol to generate a consensus 16S rRNA gene sequence from forward and reverse sequences.

DOI

dx.doi.org/10.17504/protocols.io.6qpvr49e3gmk/v1

PROTOCOL CITATION

Kenneth Acosta 2022. 16S_PCR_Sequence_Analysis. **protocols.io** https://dx.doi.org/10.17504/protocols.io.6qpvr49e3gmk/v1

KEYWORDS

Genotyping

LICENSE

This is an open access protocol distributed under the terms of the <u>Creative Commons</u>

<u>Attribution License</u>, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 09, 2022

LAST MODIFIED

Dec 12, 2022

PROTOCOL INTEGER ID

73758



1

Citation: Kenneth Acosta 16S_PCR_Sequence_Analysis https://dx.doi.org/10.17504/protocols.io.6qpvr49e3gmk/v1

Software

Serial Cloner NAME

Franck Perez DEVELOPER

Software

FinchTV NAME

macOS

Geospiza Research Team

DEVELOPER

0.S

Adding U515 sequence feature

1 Open Serial Cloner.

Software

Serial Cloner

Franck Perez

- 2 Go to "Features" → "Manage Features".
- 3 Add new collection. Name "16S rRNA Gene Conserved Regions".
- 4 Add new feature. Name "U515". Add sequence "GTGCCAGCAGCCGCGGTAA".

protocols.io

2

Citation: Kenneth Acosta 16S_PCR_Sequence_Analysis https://dx.doi.org/10.17504/protocols.io.6qpvr49e3gmk/v1

5 Check "Scan For This Feature" box. Click Ok.

	Generating consensus 16S rNA gene sequence
6	Upload 16S rRNA gene forward and reverse sequences for the sample.
7	Go to "Function" → "Align Two Sequences".
8	Check "anti-parallel" box for reverse sequence.
9	Click "Local Align".
10	Click "Build Consensus". New dialog box will appear.
11	Select "Unsaved Consensus" window then click "Features" → "Scan Sequence". "U515" feature should be highlighted in consensus sequence.
12	Click "Features" tab in the "Unsaved Consensus" window.
13	The nucleotide positions of the "U515" marker will be given. Select the "from" box in top right corner of window

Citation: Kenneth Acosta 16S_PCR_Sequence_Analysis https://dx.doi.org/10.17504/protocols.io.6qpvr49e3gmk/v1

motocols.io

3

and type in the nucleotide position that is -385 nt from the beginning of the "U515" marker. Select the "to" box in top right corner of window and type in the nucleotide position that is +216 nt from the end of the "U515" marker. A 620 bp consensus sequence should be generated.

To check for quality of consensus sequence select "Align Two Sequences" window and look for differences in base calls between forward and reverse sequences within the consensus sequence. Sometimes there will be a mismatch or gap between forward and reverse sequences. To assist in making the right base call at a certain nucleotide position use FinchTV to view chromatograms and determine the right call using the higher-quality base call.

Software	
FinchTV	NAME
macOS	OS
Geospiza Research Team	DEVELOPER

- 15 Save consensus 16S rRNA gene sequence.
- Annotate consensus 16S rRNA gene sequence by searching a database (e.g. NCBI) or using a 16S rRNA gene classifier (e.g. RDP classifier).