



Oct 18, 2021

# 🌐 Preparation of BG-11 media Agar

**Ashwinuday**<sup>1</sup><sup>1</sup>Indian Institute of Science Education and Research, Pune[dx.doi.org/10.17504/protocols.io.by6rpzd6](https://dx.doi.org/10.17504/protocols.io.by6rpzd6)**iGEM IISER Pune India 2021****Ashwinuday**

This is for making solid media of BG-11

DOI

[dx.doi.org/10.17504/protocols.io.by6rpzd6](https://dx.doi.org/10.17504/protocols.io.by6rpzd6)Ashwinuday 2021. Preparation of BG-11 media Agar . **protocols.io**<https://dx.doi.org/10.17504/protocols.io.by6rpzd6>

\_\_\_\_\_ protocol ,

Oct 18, 2021

Oct 18, 2021

54193

- 1 To make BG-11 Agar plates, add Agar-Agar Type-1 equivalent to 1.5% (w/v) of liquid media and autoclave.
- 2 Then before it solidifies (when the flask is bearable to touch), add the required antibiotics to the media.
- 3 Pour ~30-40 ml into plates soon after.

- 4 Keep the plates opened in the hood.
- 5 Once it solidifies in the plates then close those and keep the plates inverted into a cold storage.