

Sep 02, 2021

3.4 Genome Editing with CRISPR/Cas9

✓ Book Chapter

📁 In 1 collection

Evelyn J. Sauter¹, Lisa K. Kutsche¹, Simon D. Klapper¹, Volker Busskamp¹¹Center for Regenerative Therapies, Technische Universität Dresden, Dresden, Germany

1 Works for me

🔗 Share

dx.doi.org/10.17504/protocols.io.bqhemt3e

Springer Nature Books

satyavati Kharde

ABSTRACT

This is part 3.4 of the "[Induced Neurons for the Study of Neurodegenerative and Neurodevelopmental Disorders](#)" collection of protocols.

Collection Abstract: Patient-derived or genomically modified human induced pluripotent stem cells (iPSCs) offer the opportunity to study neurodevelopmental and neurodegenerative disorders. Overexpression of certain neurogenic transcription factors (TFs) in iPSCs can induce efficient differentiation into homogeneous populations of the disease-relevant neuronal cell types. Here we provide protocols for genomic manipulations of iPSCs by CRISPR/Cas9. We also introduce two methods, based on lentiviral delivery and the piggyBac transposon system, to stably integrate neurogenic TFs into human iPSCs. Furthermore, we describe the TF-mediated neuronal differentiation and maturation in combination with astrocyte cocultures.

ATTACHMENTS

[Sauter2019_Protocol_InducedNeuronsForTheStudyOfNeu.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.bqhemt3e

EXTERNAL LINK

https://link.springer.com/protocol/10.1007/978-1-4939-9080-1_9

PROTOCOL CITATION

Evelyn J. Sauter, Lisa K. Kutsche, Simon D. Klapper, Volker Busskamp 2021. 3.4 Genome Editing with CRISPR/Cas9. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bqhemt3e>

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Sauter E.J., Kutsche L.K., Klapper S.D., Busskamp V. (2019) Induced Neurons for the Study of Neurodegenerative and Neurodevelopmental Disorders. In: Ben-Yosef D., Mayshar Y. (eds) Fragile-X Syndrome. Methods in Molecular Biology, vol 1942. Humana Press, New York, NY.
https://doi.org/10.1007/978-1-4939-9080-1_9


COLLECTIONS ⓘ

**Induced Neurons for the Study of Neurodegenerative and Neurodevelopmental Disorders**

KEYWORDS

Human induced pluripotent stem cells, Nucleofection, PiggyBac transposon, Lentiviral transduction, CRISPR/Cas9, Transcription factor-mediated neuronal differentiation, Astrocyte coculture

LICENSE

 This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 07, 2020

LAST MODIFIED

Sep 02, 2021

OWNERSHIP HISTORY

Dec 07, 2020  Lenny Teytelman protocols.io

Jul 05, 2021  Emma Ganley protocols.io

Aug 24, 2021  Satyavati Kharde

Aug 26, 2021  satyavati Kharde

PROTOCOL INTEGER ID

45318

PARENT PROTOCOLS

Part of collection

[Induced Neurons for the Study of Neurodegenerative and Neurodevelopmental Disorders](#)

2. Materials

2.4 Genome Editing with CRISPR/Cas9

1. Cas9-sgRNA construct, such as pSpCas9(BB)-2A-Puro (PX459) V2.0 (Addgene plasmid #62988) or pSpCas9(BB)-2A-GFP (PX458) (Addgene plasmid #48138), which is expressing both sgRNA together with the gRNA scaffold as well as the staphylococcus pyogenes Cas9 [10]. Store at -20°C .
2. **10 x T4 Ligation Buffer (New England Biolabs)**. Store at -20°C .
3. T4 Polynucleotide Kinase (New England Biolabs). Store at -20°C .
4. BbsI restriction enzyme (**10 U/ μl**) (Thermo Fisher Scientific). Store at -20°C .
5. Buffer G (Thermo Fisher Scientific). Store at -20°C .
6. Calf Intestinal Alkaline Phosphatase (CIP) (New England Biolabs). Store at -20°C .
7. Electrophoresis gel and chamber.
8. Gel extraction kit, such as QIAquick® Gel Extraction Kit (Qiagen).
9. Ligation kit, such as Mighty Mix ligation kit (Clontech). Store at -20°C .
10. Chemically competent bacteria (e.g., Stbl3). Store at -80°C .
11. S.O.C. medium. Store at 4°C .
12. LB-Antibiotics plates with **100 $\mu\text{g}/\text{mL}$ ampicillin**. Store at 4°C .
13. Miniprep kit, such as QIAprep® Spin Miniprep Kit (Qiagen).
14. Sequencing primer, such as 5'-TTTCTTGGGTAGTTTGCAGTTT-3'. Dilute in ddH₂O to a concentration of **10 Micromolar (μM)** and store at -20°C .
15. 4D-Nucleofector™ System (see Materials section in [protocol 3.2](#)).

SAFETY WARNINGS

For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

BEFORE STARTING



NB Introduction, Notes, and References are in the [Collection Guidelines tab](#)


3.4 Genome Editing with CRISPR/Cas9 1d 7h 52m 45s

- 1 Design a sgRNA for the locus of interest using designing tools (such as <http://crispr.mit.edu/>, see **Note 23**) [12].
- 2 Order two oligos representing the sgRNA with a design as follows: top oligo—CACC(G)[20 N of sgRNA], bottom oligo—AAAC[20 N reverse complement of sgRNA](C) [10] (see **Note 24**).

3 

35m

Anneal the oligos in a thermocycler using **2 μl top oligo**, **2 μl bottom oligo**, **2 μl 10x T4 Ligation Buffer**, T4 Polynucleotide Kinase, and **13 μl ddH₂O** and program the thermocycler with  **00:30:00** at 37°C ,  **00:05:00** at 95°C , a ramp down to 25°C (ramp rate of $5^{\circ}\text{C}/\text{min}$) and hold at 4°C . The oligo hybrid can directly be used for cloning into the linearized PX459 or PX458 vector.

4 

Cut and dephosphorylate the vector using **2 µl BbsI** with **3 µg vector** and fill up to 50 µl/reaction volume with ddH₂O (minimum of 0.5 µl BbsI for 5 µg of vector).

1h 15m



Incubate for **01:00:00** at **37 °C** and optionally heat-inactivate the digestion for **00:15:00** at **65 °C**.

- 6 Gel-purify the vector by running on an agarose gel and extracting the band using a gel extraction kit, such as QIAquick[®] Gel Extraction Kit (Qiagen), according to the manufacturer's instructions. The vector can be stored at **-20 °C** for several months.



30m

Ligate sgRNA insert and vector in a 10 µl reaction with **100 ng vector**, **2 µl oligo hybrid (1:250 diluted)**, **5 µl Mighty Mix**, and ddH₂O. Incubate the mixture for **00:30:00** at **16 °C**.

- 8 Transform **5 µl reaction product** into chemically competent bacteria, such as Stbl3. The rest of the reaction can be stored for up to 2 weeks.



32m 45s

Thaw the bacteria for ~ **00:10:00** **On ice**, add the ligation reaction, stir carefully and incubate **On ice** for **00:05:00** – **00:10:00**. Perform a heat shock for **00:00:45** at **42 °C** transfer the sample back to ice for **00:02:00** – **00:05:00**. Add **300 µl S.O.C. medium** and shake the bacteria at **300 rpm, 37°C, 01:00:00**.



1h

Plate the bacteria on LB-plates with **100 µg/mL ampicillin** and incubate **Overnight** at **37 °C** and proceed with the plasmid preparation using a Miniprep kit such as the QIAprep[®] Spin Miniprep Kit (Qiagen) according to the manufacturer's guidelines.

- 11 Check the construct for correctness using the respective sequencing primer and expand the DNA to a high concentration stock.

Optionally: Test the cutting efficiency with a T7 endonuclease assay [13] (*see* **Note 25**).

- 12 If homologous recombination of a disease correction or knockin of a reporter is the aim of the gene editing, a donor construct has to be provided in addition to the Cas9-sgRNA vector(s). The tag or the gene for the knockin should be framed by regions of the locus in approximately the same length as the part that is transferred into the gene (such as ~1000 bp for incorporation of an antibiotic resistance with promoter) and provided as a linearized DNA fragment. Alternatively, 5–25 bp homology can be used with microhomologymediated end-joining [14].

- 13 Electroporate **10 µg vector** encoding for the correct sgRNA and optionally ~ **5 µg donor construct** into iPSCs (*see* materials and steps in [protocol 3.2](#)). Seed the cells into two wells of a 6-well plate.

- 14 Using the PX459 vector, it is advisable to use [M]0.5 µg/mL – [M]0.8 µg/mL puromycin for ⌚ 12:00:00 – 1d 4h ⌚ 16:00:00 for the next night after electroporation in combination with ROCKi.
- 15 Keep ROCKi for 2 days and check if you obtained single cells that grow to small colonies. Once the colonies are visible by eye, check half of a colony for integration or knockout of the target gene by (colony) PCR with primers binding around or in the locus that you tested for specificity before (using control DNA). Transfer the other part of each colony not used for PCR to a new well (of a 48-well plate). These are the monoclonal knockout lines usable if the PCR shows the respective shifts.
- 16 Once a potentially positive clone has been detected, sequence the locus using Sanger sequencing, optimally with subcloning of the PCR product into a carrier vector (for instance using TOPO cloning).