



Mar 01, 2022

🌐 Quick Decapitation for Paraffin Embedding

Haley Geertsma¹¹University of Ottawa

1

dx.doi.org/10.17504/protocols.io.b5szq6f6

Haley Geertsma
University of Ottawa

This protocol is used to quickly decapitate mice and isolate their brain for paraffin embedding and microtome sectioning.

DOI

dx.doi.org/10.17504/protocols.io.b5szq6f6

Haley Geertsma 2022. Quick Decapitation for Paraffin Embedding. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b5szq6f6>



_____ protocol ,

Mar 01, 2022

Mar 01, 2022

58937

- 1 Sedate mice with isoflurane then quickly decapitate. 5m
- 2 Extract the brain and incubate in 4% paraformaldehyde or 10% formalin for 72 hours 3d
- 3 Incubate brain in 70% ethanol and submit for paraffin embedding. Once embedded, section^{1h} brains at 5µm.

