


Aug 19, 2022

ACE-inhibitory activity assay: IC₅₀

Ezequiel Coscueta¹, & Maria Manuela Pintado¹¹Centro de Biotecnologia e Química Fina (CBQF)-Laboratório Associado, Escola Superior de Biotecnologia, Universidade Católica Portuguesa, Rua Diogo Botelho 1327, 4169-005 Porto, Portuga

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.q26g74q5kgwz/v1

Ezequiel Coscueta

ABSTRACT

This protocol describes the procedure for the determination of the IC₅₀ in inhibition on angiotensin-I converting enzyme (ACE) activity. ACE is also known as peptidyl dipeptidase A because it removes C-terminal dipeptides from a wide variety of peptide substrates. In the assay described here, the chosen substrate is the intramolecularly quenched fluorescent tripeptide o-aminobenzoylglycyl-p-nitrophenylalanylproline (Abz–Gly–Phe(NO₂)–Pro). Hydrolysis of this substrate by the action of ACE generates the fluorescent product o-aminobenzoylglycine (Abz–Gly).

DOI

dx.doi.org/10.17504/protocols.io.q26g74q5kgwz/v1

PROTOCOL CITATION

Ezequiel Coscueta, & Maria Manuela Pintado 2022. ACE-inhibitory activity assay: IC₅₀. **protocols.io**
<https://protocols.io/view/ace-inhibitory-activity-assay-ic50-cfejtjcn>



MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Coscueta, E. R., Campos, D. A., Osório, H., Nerli, B. B., & Pintado, M. (2019). Enzymatic soy protein hydrolysis: A tool for biofunctional food ingredient production. Food Chemistry: X, 1(100006). <https://doi.org/10.1016/j.fochx.2019.100006>

KEYWORDS

ACE, angiotensin-converting enzyme, antihypertensive, IC₅₀, bioactivity, ngiotensin-I converting enzyme

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 17, 2022

LAST MODIFIED

Aug 19, 2022

PROTOCOL INTEGER ID

68779

GUIDELINES

ACE-inhibitory activity is measured by fluorescence using the method of Sentandreu & Toldrá (2006) with some modifications (Coscueta, Campos, Osório, Nerli, & Pintado, 2019). A total of 40 µL of ultrapure water or ACE working solution are added to each microtiter-plate well, then adjusted to 80 µL by adding ultrapure water to blank, control or samples. For direct or 1/2 diluted samples a sample blank is also made. The enzyme reaction is started by adding 160 µL of substrate solution and the mixture is incubated at 37 °C. Serial dilutions are made of each sample, usually 1/1 to 1/32. The fluorescence generated is measured at 30 min using a Multidetector plate reader (Synergy H1, Vermont, USA). The assay is performed in a black 96-well microplate (Nunc, Denmark). Excitation and emission wavelengths are 350 and 420 nm, respectively.

Note: In some cases it is necessary to extend the range to 1/5 or 1/10 serial dilutions.

Note: Generally the most concentrated dilutions are the ones that can fluoresce as the sample itself, so in the example sample blanks (SPLB) are used only for the two most concentrated dilutions. In the case of highly fluorescent samples, SPLB should be used for all dilutions.

Sentandreu, M. Á., & Toldrá, F. (2006). A rapid, simple and sensitive fluorescence method for the assay of angiotensin-I converting enzyme. Food Chemistry.
<https://doi.org/10.1016/j.foodchem.2005.06.006>

Coscueta, E. R., Campos, D. A., Osório, H., Nerli, B. B., & Pintado, M. (2019). Enzymatic soy protein hydrolysis: A tool for biofunctional food ingredient production. Food Chemistry: X.
<https://doi.org/10.1016/j.fochx.2019.100006>

MATERIALS TEXT

[Angiotensin Converting Enzyme from rabbit lung](#) **Sigma** –

Aldrich Catalog #A6778

[Tris Base](#) **Fisher**

Scientific Catalog #BP152-1

or [Tris-HCl](#) **Sigma Aldrich**

[Ultrapure Water](#) **Contributed by users**

[ZnCl₂ anhydrous](#) **Thermo Fisher**

Scientific Catalog #11497737

[Sodium Chloride](#) **Fisher**

Scientific Catalog #S271

[Abz-Gly-p-nitro-Phe-Pro-OH](#)

[Trifluoroacetate](#) **bachem Catalog #4003531.0050**

Synergy H1M

Multidetector plate reader

BioTek 679SH1M-SN

Detection modes: UV-Vis absorbance;
fluorescence intensity; luminescence

Read methods: endpoint, kinetic, spectral
scanning, well area scanning

Microplate types: 6- to 384-well plates

Temperature control: to 45 °C

Shaking: linear, orbital, double orbital

96 Well Black Plate

microplate

Nunc 237108

96-Well black polystyrene microplate, flat
bottom, 50 to 250 µL, non-sterile, non-treated

96 Well Black Plate

BEFORE STARTING

Prepare the necessary reagents carefully.

[M]1 U/mL ACE stock solution

- Dissolve the

⌘ [Angiotensin Converting Enzyme from rabbit lung Sigma –](#)

[Aldrich Catalog #A6778](#)

(peptidyl-dipeptidase A, EC 3.4.15.1) in a solution of 50% glycerol in ultrapure water, to obtain a final concentration of 1 U/mL.

- Make aliquots of **200 µL** of the solution and store at **-20 °C**.

[M]0.150 Molarity (M) Tris buffer **pH8.3**

- Dissolve **1.817 g** of Tris base (MW = 121.14) or **2.364 g** Tris-HCl (MW = 157.60) in approx. **90 mL** ultrapure water.
 - Titrate to **pH8.77** at the lab temperature of **Room temperature** with monovalent strong base or acid as needed.
 - Make up volume to **100 mL** with ultrapure water.
- Buffer will be **pH8.3** at **37 °C**.

[M]42 mU/mL ACE work solution

- Prepare a solution of **1 millimolar (mM) ZnCl₂**, dissolving **1.4 mg** ZnCl₂ (MW = 136.28) in **10 mL** ultrapure water. Store at **-20 °C**.
- Prepare **0.1 millimolar (mM) ZnCl₂** **0.150 Molarity (M) Tris buffer** **pH8.3** (Enzyme buffer). Dilute 1/10 the previous solution (**0.1 millimolar (mM)**) and add **25 µL** of this solution to **25 mL** of **0.150 Molarity (M) Tris buffer** **pH8.3**. Store at **4 °C** for a maximum of one week, or at **-20 °C** for a maximum of six months.
- Dilute 1/24 the **1 U/mL ACE stock solution** with the Enzyme buffer. Prepare daily.

[M]0.45 millimolar (mM) Substrate solution

- Prepare **1.125 Molarity (M) NaCl** **0.150 Molarity (M) Tris buffer** **pH8.3** (Substrate buffer). Dissolve **3.2872 g** NaCl (MW = 58.44) in **50 mL** of **0.150 Molarity (M) Tris buffer** **pH8.3**. Store at **4 °C** for a maximum of one week, or at **-20 °C** for a maximum of six months.
- Dissolve 3.6 mg of substrate [⌘ Abz-Gly-p-nitro-Phe-Pro-OH](#) [Trifluoroacetate](#) [bachem Catalog #4003531.0050](#) in **16 mL** Substrate buffer (for 96 wells). Prepare at the moment.

Analysis

30m

1

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	BLK	BLK	SPLB1:1 1	SPLB1:1 1	SPLB2:1 1	SPLB2:1 1	SPLB3:1 1	SPLB3:1 1	SPLB4:1 1	SPLB4:1 1	
B	CTL	CTL	CTL	SPLB1:2 2	SPLB1:2 2	SPLB2:2 2	SPLB2:2 2	SPLB3:2 2	SPLB3:2 2	SPLB4:2 2	SPLB4:2 2	
C	SPL1:1 1	SPL1:1 1	SPL1:1 1	SPL2:1 1	SPL2:1 1	SPL2:1 1	SPL3:1 1	SPL3:1 1	SPL3:1 1	SPL4:1 1	SPL4:1 1	SPL4:1 1
D	SPL1:2 2	SPL1:2 2	SPL1:2 2	SPL2:2 2	SPL2:2 2	SPL2:2 2	SPL3:2 2	SPL3:2 2	SPL3:2 2	SPL4:2 2	SPL4:2 2	SPL4:2 2
E	SPL1:3 4	SPL1:3 4	SPL1:3 4	SPL2:3 4	SPL2:3 4	SPL2:3 4	SPL3:3 4	SPL3:3 4	SPL3:3 4	SPL4:3 4	SPL4:3 4	SPL4:3 4
F	SPL1:4 8	SPL1:4 8	SPL1:4 8	SPL2:4 8	SPL2:4 8	SPL2:4 8	SPL3:4 8	SPL3:4 8	SPL3:4 8	SPL4:4 8	SPL4:4 8	SPL4:4 8
G	SPL1:5 16	SPL1:5 16	SPL1:5 16	SPL2:5 16	SPL2:5 16	SPL2:5 16	SPL3:5 16	SPL3:5 16	SPL3:5 16	SPL4:5 16	SPL4:5 16	SPL4:5 16
H	SPL1:6 32	SPL1:6 32	SPL1:6 32	SPL2:6 32	SPL2:6 32	SPL2:6 32	SPL3:6 32	SPL3:6 32	SPL3:6 32	SPL4:6 32	SPL4:6 32	SPL4:6 32

BLK Blank

CTL Control

SPLB Sample blank

SPL Sample

Microplate outline example for IC₅₀ determination.



1.1 Add  **80 µL ultrapure water** to blank (BLK)

1.2 Add  **40 µL ultrapure water** to control (CTL) and sample blank (SPLB)

1.3 Add  **40 µL sample dilution** to sample (SPL) and sample blank (SPLB)

1.4 Add  **40 µL ACE work solution** to control (CTL) and sample (SPL)

1.5 Add  **160 µL Substrate solution** to control all wells


1.6 Incubate at  **37 °C**  **00:30:00** and read fluorescence with **350 nm excitation**^{30m} wavelength and **420 nm emission** wavelength

Result treatment

- 2 Inhibitory activity is expressed as the peptide concentration required to inhibit the original ACE activity by 50% (IC_{50}). The formula applied to calculate de percentage of ACE-inhibitory is:

$$iACE (\%) = ((F_{CTL} - F_{BLK}) - (F_{SPL} - F_{SPLB})) * \frac{100}{F_{CTL} - F_{BLK}}$$

Non-linear fitting to the data is performed to calculate the IC_{50} values, using the 5 Parameter curve fit method and then Interpolating to 50.



Example of a typical inhibition curve as a function of dilution factor.