



FEB 28, 2024

🌐 Odor Retention

sdwalto¹, Jeffrey Kordower¹, Bryan_Killinger²

¹ASU NDRC; ²Rush University



sdwalto

ABSTRACT

Odor retention test optimized for mice. This test assesses short term olfactory memory. Mice with olfactory impairments or deficits should score lower.

GUIDELINES

ODERANT PAIRS: limonene (diluted 1:10) and carvone (diluted 1:10); Amylacetate (1:100) and Anisol (1:100); and propyl acetate (1:10) and Pentanal (1:10).

Notes:

- Paper swabs soaked in a given odorant were used inside the cartridges.
- Investigation time defined as the duration of active sniffing with the nose placed less than 1 cm away from the cartridge.

MATERIALS

- ODERANT PAIRS: limonene (diluted 1:10) and carvone (diluted 1:10); Amylacetate (1:100) and Anisol (1:100); and propyl acetate (1:10) and Pentanal (1:10).

- Paper swabs

- Mesh balls (tea balls)

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.eq2lyjoeqlx9/v1

Protocol Citation: sdwalto, Jeffrey Kordower, Bryan_Killinger 2024. Odor Retention.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.eq2lyjoeqlx9/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Feb 28, 2024

Last Modified: Feb 28, 2024

PROTOCOL integer ID: 95895

Keywords: ASAPCRN, odor threshold, mice behavior

Funders Acknowledgement:

Jeffrey Kordower

Grant ID: NIH R21 NS109871

- 1 Prepare and empty rodent cage with no bedding.
- 2 Prepare the odorants by diluting in mineral oil. (See guidelines)
- 3 Acclimate mouse to setup (empty cage with odor cartridges). Place mouse in the cage and expose to 2 empty odor cartridges on either side of the cage for 5 mins.
- 4 Expose mouse to the same unfamiliar odor in both cartridges for 5 mins. (**Acquisition**)
- 5 After 6, 16 or 30 mins post-acquisition, expose the mouse to a novel odorant and the familiar odorant from Step 4 and record the time spent investigating both the novel and familiar odorant for all trials. (Total investigation time capped at 1 min). (**Recall**)
- 6 The preferred odor was recorded and the preference was calculated based on the recall trial.

