

Jul 10, 2024

Primary cortical astrocyte isolation and culture

DOI

dx.doi.org/10.17504/protocols.io.x54v9dqqmg3e/v1

Shiyi Wang¹

¹Duke University

ASAP Collaborative Rese...



Shiyi Wang

Duke University

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.x54v9dqqmg3e/v1**

Document Citation: Shiyi Wang 2024. Primary cortical astrocyte isolation and culture. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.x54v9dqqmg3e/v1>

License: This is an open access document distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Created: May 23, 2023

Last Modified: July 10, 2024

Document Integer ID: 82322

Keywords: ASAPCRN

Funders Acknowledgement:
Aligning Science Across
Parkinson's (ASAP) initiative
Grant ID: ASAP-020607

Abstract

Primary cortical astrocyte isolation and culture



1. P1 rat cortices from both sexes were micro-dissected, papain digested, triturated in low and high ovomucoid solutions, and resuspended in astrocyte growth media (AGM: DMEM (GIBCO 11960), 10% FBS, 10 μ M, hydrocortisone, 100 U/ml Pen/Strep, 2 mM L-Glutamine, 5 μ g/ml Insulin, 1 mM Na Pyruvate, 5 μ g/ml N-Acetyl-L-cysteine).
2. Between 15-20 million cells were plated on 75 mm² flasks (non-ventilated cap) coated with poly-D-lysine and incubated at 37°C in 10% CO₂.
3. On DIV 3, the removal of non-astrocyte cells was performed by forcefully shaking closed flasks by hand for 10-15 s until only an adherent monolayer of astrocytes remained.
4. AraC was added to the media from DIV 5 to DIV 7 to eliminate contaminating fibroblasts. On DIV 7, astrocytes were trypsinized (0.05% Trypsin-EDTA) and plated into 12-well or 6-well dishes.
5. On DIV 8, cultured rat astrocytes were transfected with shRNA and/or expression plasmids using Lipofectamine LTX with Plus Reagent (Thermo Scientific) per the manufacturer's protocol.
6. Briefly, 1 μ g (12-well) or 2 μ g (6-well) total DNA was diluted in Opti-MEM containing Plus Reagent, mixed with Opti-MEM containing LTX (1:2 DNA to LTX), and incubated for 30 minutes.
7. The transfection solution was added to astrocyte cultures and incubated at 37°C for 3 hours. On DIV 10, astrocytes were trypsinized, resuspended in NGM plus, plated (20,000 cells per well) onto DIV 10 neurons, and co-cultured for 48 hours.