

Sep 09, 2024

Fluorescent Gelatin Degradation Assay to Evaluate EVh Action in TME Cells

This protocol is a draft, published without a DOI.

Bianca Cruz Pachane¹, Heloisa Sobreiro Selistre de Araujo¹

¹Department of Physiological Sciences, Universidade Federal de São Carlos - UFSCar



Bianca Cruz Pachane

Universidade Federal de São Carlos - UFSCar

OPEN  ACCESS



Protocol Citation: Bianca Cruz Pachane, Heloisa Sobreiro Selistre de Araujo 2024. Fluorescent Gelatin Degradation Assay to Evaluate EVh Action in TME Cells. **protocols.io** <https://protocols.io/view/fluorescent-gelatin-degradation-assay-to-evaluate-dkt34wqn>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 07, 2024

Last Modified: September 09, 2024

Protocol Integer ID: 107099

Keywords: Cell Invasion, Gelatin Degradation, Gelatinase

Funders Acknowledgement:

**São Paulo Research
Foundation**

Grant ID: 2021/01983-4



Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

The fluorescent gelatin degradation assay is a method to study cell invasion by detecting gelatinase activity *in vitro* upon epifluorescence microscopy analysis. In this protocol, the method has been applied to evaluate the effect of hypoxic EVs from TNBC cell line MDA-MB-231 in four cellular models for the tumor microenvironment - MDA-MB-231 (tumor cell), HUVEC (endothelial cell), HDFa (dermal fibroblast) and THP-1 (monocyte). Adapted from Pachane et al (2022) (PMID: 36293503).

Materials

Materials and reagents

1. Corning 96-well Flat Clear Bottom Black Polystyrene TC-treated Microplates, Individually Wrapped
2. Sterile microtubes and pipettes
3. Gelatin From Pig Skin, Fluorescein Conjugate, Thermo Fisher - Catalog #G13187
4. Sterile PBS
5. OptiMEM I Reduced Serum Media, Gibco - Catalog #31985070
6. Trypan Blue solution 0.4%, Merck Millipore (Sigma-Aldrich) - Catalog #T8154-100 ml
7. Paraformaldehyde solution (PFA 4% in deionized water, pH 7.6 - Sterile)
8. Triton X-100 0.1% (v/v) in deionized water
9. Phalloidin + DAPI (1 μ l Phalloidin-iFluor 647, Abcam - Catalog #ab176759 + 0.76 μ l 4,6-Diamidino-2-Phenylindole, Dihydrochloride (DAPI), Thermo Fisher Scientific - Catalog #D1306 in 5 mL PBS)

Cell lines and growth media

- MDA-MB-231 (ATCC® CRM-HTB-26™) - Leibovitz L-15 10% FBS
- HDFa (ATCC® PCS-201-012™) - DMEM 10% FBS 1% pen/strep
- HUVEC (ATCC® CRL-1730™) - DMEM 10% FBS 1% pen/strep
- THP-1 (ATCC® TIB-202™) - RPMI 1640 10% FBS 1% pen/strep

Equipments:

1. Biological cabinet
2. Cell incubator (37 °C, 5% CO₂)
3. Cell counter - TC20 Cell Counter, Bio-Rad - Catalog #1450011
4. Epifluorescence microscope - ImageXpress Micro XLS, Molecular Devices - Catalog #500496



Protocol materials



Parafilm

Step 18

Safety warnings



Light-sensitive assay. Work under sterile conditions.

Before start

Fluorescent gelatin preparation: Under sterile conditions, solubilize the fluorescent gelatin stock at $37\text{ }^{\circ}\text{C}$ with warmed PBS following the manufacturer's instructions for a concentration of 5 mg/mL . Aliquot in microtubes and maintain at $-20\text{ }^{\circ}\text{C}$ until time of use.

Before use, thaw gelatin at $37\text{ }^{\circ}\text{C}$ for $00:30:00$. Dilute stock to a 0.2 mg/mL working solution with warmed PBS and maintain at $37\text{ }^{\circ}\text{C}$ until use.

Cell culture: Maintain cells in culture during at least two passages after thawing.

Fluorescent Gelatin Coating

30m

- 1 Open a new 96-well black plate under sterile conditions and label groups in technical triplicates to contain a **vehicle (PBS) control** (i.e., untreated cells in OptiMEM) and the **EVh-treated** group (i.e., EVh-treated cells in OptiMEM) for each cell line.

1.1 Experimental plate map:

	1	2	3	4	5	
A	MDA-MB-231 P BS	MDA-MB-231 P BS	MDA-MB-231 P BS	MDA-MB-231 E Vh	MDA-MB-231 E Vh	MDA-MB-231 E Vh
B	HUVEC PBS	HUVEC PBS	HUVEC PBS	HUVEC EVh	HUVEC EVh	HUVEC EVh
C	HDFa PBS	HDFa PBS	HDFa PBS	HDFa EVh	HDFa EVh	HDFa EVh
D	THP-1 PBS	THP-1 PBS	THP-1 PBS	THP-1 EVh	THP-1 EVh	THP-1 EVh
E						
F						
G						
H						



	7	8	9	10	11	
A						
B						
C						
D						
E						
F						
G						
H						

2 Apply 70 μ L of the fluorescent gelatin working solution [M] 0.2 mg/mL directly to the bottom of each well and prevent the formation of bubbles.

3 Incubate plate for 00:30:00 at 37 °C 5% CO₂ .

30m





4 Carefully remove excess coating (avoid touching well bottom).

5 Pre-condition coating with  200 μ L OptiMEM for  00:30:00 at  37 °C 5% CO₂ .

30m



Cell seeding

1d


6 Subculture cells as usual. Resuspend cell pellet in OptiMEM and count cells using the trypan blue exclusion method.

7 Remove pre-conditioning media from the wells (avoid touching well bottom).

8 Add cell suspension into each well to a total volume of 200 μ L:

- MDA-MB-231: 5×10^3 cells/well (= 1×10^5 células/ml)
- HUVEC: 5×10^3 cells/well (= 1×10^5 células/ml)
- HDFa: 2×10^3 cells/well (= 1×10^4 células/ml)
- THP-1: 5×10^3 cells/well (= 1×10^5 células/ml)

9 Treat cells with EVh (10^9 particles/ml) or the equivalent treatment volume in PBS.

10 Incubate plate for  24:00:00 at  37 °C 5% CO₂

1d




Fixation and Cell Staining

10m




11 Remove the supernatant by aspiration.

12 Fix cells with  100 μ L warmed 4% PFA at  Room temperature for  00:10:00

10m









13 Wash wells twice with  100 μ L PBS

II







14 Permeabilize cells with  100 μ L 0.1% Triton X-100 at  Room temperature for  00:05:00

5m





- 15 Wash wells twice with  100 μ L PBS
- 16 Stain cells with the DAPI + Phalloidin-647 mixture. Add  100 μ L of staining solution to each well and incubate at  Room temperature , protected from light for  00:20:00 . 20m
- 17 Wash wells twice with  100 μ L PBS
- 18 Maintain wells with  200 μ L PBS for analysis. Seal the plate with  Parafilm **Contributed by users** and cover it with aluminum foil for storage at  4 $^{\circ}$ C for up to 6 months. II

Cell Imaging by Epifluorescence HTS

- 19 Using the microscope **ImageXpress Micro XLS+ (Molecular Devices)**, check the template for the Corning 3603 plate and the filters for DAPI (nuclei), FITC (gelatin) and Cy5 (phalloidin-647). 
- 20 Set laser intensity to a minimum of 10 ms and increase gradatively if necessary. 
- 21 Check the wells using the 4X objective. 
- 22 Change into the 20x objective and adjust the laser focus. Select 9 sites per well minimally. 
- 23 Acquire the plate. Export metadata for analysis. 
- 24 For representative images, change into the 40x objective and adjust the laser focus. Select the sites of interest and acquire. Export image channels and combinations. 


Gelatin Degradation Quantification on FIJI

- 25 On FIJI (ImageJ), import HTD files through BioFormats. 
- 26 Images should already be scaled. If not, adjust scale based on the objective lens used for acquisition. 




- 27 Set measurements to contain "Area", "Standard Deviation", "Shape Descriptor", "Mean grey value", "Perimeter" and "Display label".
- 28 Concatenate all stacks into a single hyperstack.
- 29 Split channels and select the FITC stack for analysis.
- 30 Set a threshold to encompass the degradation spots but not the background. Write down the threshold values. Create a new stack with the binary images.
- 31 To measure the degraded area, analyze particles with a range of "5-Infinity" and select "Summarize".
- 32 Save CSV file. The degraded area (in μm^2) per site will be compared between groups in the statistical analysis.

Cell Counting on FIJI

- 33  [go to step #25](#) and follow through step #27
- 34 Split channels and select the DAPI stack for counting.
- 35 Set a threshold to contain nuclei. Create a new stack with the binary images.
- 36 To count cells, analyze particles with a range of "10-infinity" and select "Summarize".
- 37 Save CSV file.


Cell Morphology Analysis on FIJI

- 38  [go to step #25](#) and follow through step #27
- 39 Split channels and select the Cy5 stack for analysis.



- 40 Duplicate the stack as a guide.
- 41 Set a threshold to encompass cell cytoplasm. Create a new stack with the binary images.
- 42 Using the duplicated stack as a guide, section cells using the "pencil" tool with a 3 px thickness.
- 43 To analyze cell morphology, analyze particles with a range of "10-Infinity" and check "Clear Results".
- 44 Save CSV file. The cell circularity index of each cell will be compared between groups in the statistical analysis.

Image Processing for Representative Cells

- 45  [go to step #25](#) and follow through step #27
- 46 Split the channels of the stacks of interest.
- 47 Adjust channel colors using the "Lookup Tables" menu.
- 48 Select the cell of interest in a 200x200 px squared selection.
- 49 Save selections in each channel and the combination of all channels in PNG images.



Protocol references

PACHANE, Bianca Cruz et al. Small Extracellular Vesicles from Hypoxic Triple-Negative Breast Cancer Cells Induce Oxygen-Dependent Cell Invasion. International Journal of Molecular Sciences, [s. l.], v. 23, n. 20, p. 12646, 2022.

EVEN-RAM, Sharona; ARTYM, Vira. Extracellular Matrix Protocols: Second Edition. [S. l.]: Humana Press, 2009.