



Feb 05, 2021

Protocols for "Efficient and stable metabarcoding sequencing data using DNBSEQ-G400 sequencer validated by comprehensive community analyses"

Xiaohuan Sun¹, Yuehua Hu², Zewei Song¹¹BGI-Shenzhen; ²CAS Key Laboratory of Tropical Forest Ecology


1 Works for me dx.doi.org/10.17504/protocols.io.brtpm6mn

Xiaohuan Sun

ABSTRACT

Sequencing in DNBSEQ-G400 of this study was carried out in BGI-Shenzhen. Amplicons from 1276 soil samples were randomly and equally distributed into 3 DNBSEQ runs for both 2×200 paired-end (PE200) and 400 bases single-end sequencing (SE400). The left libraries of total 1276 samples were repetitively sequenced two more times by SE400 using different DNBSEQ-G400 machines. Three randomly chosen soil samples from each forest plot (ALS268, LJ105 and NBH217) were all amplified 3 times separately as PCR replicates. The corresponding 9 sequencing libraries were sequenced repetitively by 3 DNBSEQ runs.

Sequencing was carried out using the entire lane of the Illumina MiSeq platform in the University of Minnesota Genomic Center. Amplicons of the 1276 forest soil were normalized and pooled randomly in equal molar ratios into 8 parallel sequencing runs on MiSeq platform. Again, three randomly chosen soil samples from each forest plot (ALS268, LJ105 and NBH217) were all amplified 3 times separately as PCR replicates. The corresponding 9 sequencing libraries were sequenced repetitively by 8 MiSeq runs.



Schematic representation of sequencing strategy at DNBSEQ-G400 and Illumina MiSeq platforms.

DOI

dx.doi.org/10.17504/protocols.io.brtpm6mn

COLLECTION CITATION

Xiaohuan Sun, Yuehua Hu, Zewei Song 2021. Protocols for "Efficient and stable metabarcoding sequencing data using DNBSEQ-G400 sequencer validated by comprehensive community analyses". **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.brtpm6mn>

LICENSE

This is an open access collection distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jan 26, 2021

LAST MODIFIED

Feb 05, 2021


COLLECTION INTEGER ID

46671

ABSTRACT









Sequencing in DNBSEQ-G400 of this study was carried out in BGI-Shenzhen. Amplicons from 1276 soil samples were randomly and equally distributed into 3 DNBSEQ runs for both 2×200 paired-end (PE200) and 400 bases single-end sequencing (SE400). The left libraries of total 1276 samples were repetitively sequenced two more times by SE400 using different DNBSEQ-G400 machines. Three randomly chosen soil samples from each forest plot (ALS268, LJ105 and NBH217) were all amplified 3 times separately as PCR replicates. The corresponding 9 sequencing libraries were sequenced repetitively by 3 DNBSEQ runs.

Sequencing was carried out using the entire lane of the Illumina MiSeq platform in the University of Minnesota Genomic Center. Amplicons of the 1276 forest soil were normalized and pooled randomly in equal molar ratios into 8 parallel sequencing runs on MiSeq platform. Again, three randomly chosen soil samples from each forest plot (ALS268, LJ105 and NBH217) were all amplified 3 times separately as PCR replicates. The corresponding 9 sequencing libraries were sequenced repetitively by 8 MiSeq runs.



Schematic representation of sequencing strategy at DNBSEQ-G400 and Illumina MiSeq platforms.

FILES

-  High-throughput microbial metabarcoding amplification Protocol sequenced by DNBSEQ-G400
Version 1
by Xiaohuan Sun
-  Library construction of metabarcoding at DNBSEQ-G400 with MGIEasy universal DNA library prep Kit
Version 1
by Xiaohuan Sun
-  Sequencing at DNBSEQ-G400 platform with PE200 and SE400 modes
Version 1
by Xiaohuan Sun
-  Sample preparation for Illumina MiSeq Dual Index metabarcoding
Version 1
by Xiaohuan Sun