

MPAPASS - Creating an MPAPASS database

Joshua A Welsh¹, Sean Cook¹, Jennifer Jones¹

¹Translational Nanobiology Section, Laboratory of Pathology, Center for Cancer Research, National Cancer Institute, National Institutes of Health

Sep 19, 2020

 Works for me dx.doi.org/10.17504/protocols.io.bjnxkmfn

Translational Nanobiology Section



Joshua Welsh

ABSTRACT

This collection contains the protocols required for each step in the mpapass software pipeline for performing stitched multiplex analysis. This is one of a number of protocols in the pipeline for using the mpapass software package and is applicable to the latest release of the software.

DOI

dx.doi.org/10.17504/protocols.io.bjnxkmfn

PROTOCOL CITATION

Joshua A Welsh, Sean Cook, Jennifer Jones 2020. MPAPASS - Creating an MPAPASS database.

protocols.io

<https://dx.doi.org/10.17504/protocols.io.bjnxkmfn>



LICENSE

 This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 13, 2020

LAST MODIFIED

Sep 19, 2020

PROTOCOL INTEGER ID


40375

DISCLAIMER:

This protocol summarizes key steps for a specific type of method, which is one of a collection of methods and assays used for EV analysis in the NCI Translational Nanobiology Section at the time of submission of this protocol. Appropriate use of this protocol requires careful, cohesive integration with other methods for EV production, isolation, and characterization.

Open New Dataset Window

- 1 Open the MPAPASS software and navigate to the Menu tab in the upper left-hand corner. Under the Menu tab, choose the New Dataset option and a new window will pop-up as shown below:



Select CSV Directory

- 2 To construct the database, follow the steps outlined on the right side of the window.

First, choose a CSV Directory by clicking on the 'Select CSV Directory' button. Simply choose the folder that contains all the desired .csv files for data analysis.

The .csv files in the selected folder should now appear in the 'Directory Files' listbox on the left. Additionally, the command window should display a 'Directory Selected' message.

Bead Info


- 3 While the 'Bead Info' can be constructed as desired by using the Identifier, Marker, and Isotype textbox, the MACSPlex Human EV dataset has been automatically saved as a preset and is recommended for use with the MACSPlex Human Exosome Kit.

For constructing 'Bead Info', plug the Identifier, Marker, and Isotype into their respective boxes and then click on the 'Add Population' button. To remove a population, highlight the population and press the 'Remove Population' button. To save the new 'Bead Info' for further use, click on the 'Save Preset' button and choose a desired name.

Construct the Database

- 4 With the choice of desired 'Bead Info', it is now possible to construct the database. Using the MACSPlex Human EV preset, the screen should be similar to the one shown below.

Click on 'Construct Database' on the right and name the .xlsx file.



- 5** Outside of the MPAPASS software, navigate to the newly created .xlsx file and open it in Microsoft Excel. The file should be similar to the one shown below. Note the five worksheet tabs at the bottom of the file.

In this protocol, we will only be completing the bare minimum necessary for the database to be imported.

	A1	B	C	D	E	F	G	H	I
1	Sample_Filename_Prefix	Sample_Set_ID	Sample_ID	Sample_Grouping_ID	Sample_Control_ID	Sample_Label_Mix_No	Incubated_Sample_Volume_Microliters	Incubated_Sample_Concentration_per_ml	Sample_Source
2	export_A10_								
3	export_A11_								
4	export_A12_								
5	export_A1_								
6	export_A2_								
7	export_A3_								
8	export_A4_								
9	export_A5_								
10	export_A6_								
11	export_A7_								
12	export_A8_								
13	export_A9_								
14	export_B10_								
15	export_B11_								
16	export_B12_								
17	export_B1_								
18	export_B2_								
19	export_B3_								
20	export_B4_								
21	export_B5_								
22	export_B6_								
23	export_B7_								
24	export_B8_								
25	export_B9_								
26	export_Blank_Control_								
27	export_C10_								
28	export_C11_								
29	export_C12_								
30	export_C1_								
31	export_C2_								
32	export_C3_								
33	export_C4_								
34	export_C5_								
35	export_C6_								
36	export_C7_								
37	export_C8_								
38	export_C9_								
39	export_D10_								
40	export_D11_								
41	export_D12_								
42	export_D1_								

6 Sample Worksheet Headings:

- The first column labeled 'Sample_Filename_Prefix' gives the filename of the exported .csv files. Do not change this column unless the filename of the exported .csv file is changed.
- The second column labeled 'Sample_Set_ID' allows for grouping of files, usually by cell line and concentration.
- The third column labeled 'Sample_ID' allows for the naming of the samples. Filling this out by cell line and

concentration is encouraged as it makes data interpretation easier for the user. Any string of text is allowed. Note that the Sample ID corresponds to the Sample Set ID, thus **samples with the same Sample Set ID cannot have different Sample IDs**.

- The fourth column labeled 'Sample_Grouping_ID' allows for further classification of the sample is is related to the 'Sample_Set_ID'. Thus, in this example the Sample Set IDs of 1 and 2 refer to the Sample Grouping ID of VOK111. Similar to Sample ID, **samples with the same Sample Set ID cannot have different Sample Grouping IDs**. Although this column is recommended for understanding of the data analysis, if you do not want to use it, simply fill it out with ones.
- The fifth column labeled 'Sample_Control_ID' tells the software which file to use as a control for the sample. This is linked with the 'Controls' worksheet and will be discussed later.
- The sixth column labeled 'Sample_Label_Mix_No' tells the software what label mix was used. This is linked with the 'Labelling' worksheet and will be discussed later. **There should be no repeated Label Mix No's within a Sample Set ID as this would indicate a duplicate**.
- The other columns can be filled out as desired, but are not critical for the data analysis.

	A	B	C	D	E	F	G
	Sample_Filename_Prefix	Sample_Set_ID	Sample_ID	Sample_Grouping_ID	Sample_Control_ID	Sample_Label_Mix_No	Incubated_Sample_Volume_Microliters
1	Sample						
2	export_A10	1 VOK1111e9	VOK111	VOK111	1	10	
3	export_A11	1 VOK1111e9	VOK111	VOK111	1	11	
4	export_A12	1 VOK1111e9	VOK111	VOK111	1	12	
5	export_A1	1 VOK1111e9	VOK111	VOK111	1	1	
6	export_A2	1 VOK1111e9	VOK111	VOK111	1	2	
7	export_A3	1 VOK1111e9	VOK111	VOK111	1	3	
8	export_A4	1 VOK1111e9	VOK111	VOK111	1	4	
9	export_A5	1 VOK1111e9	VOK111	VOK111	1	5	
10	export_A6	1 VOK1111e9	VOK111	VOK111	1	6	
11	export_A7	1 VOK1111e9	VOK111	VOK111	1	7	
12	export_A8	1 VOK1111e9	VOK111	VOK111	1	8	
13	export_B0	2 VOK1111e9	VOK111	VOK111	1	9	
14	export_B10	2 VOK1111e9	VOK111	VOK111	1	10	
15	export_B11	2 VOK1111e9	VOK111	VOK111	1	11	
16	export_B12	2 VOK1111e9	VOK111	VOK111	1	12	
17	export_B1	2 VOK1111e9	VOK111	VOK111	1	1	
18	export_B2	2 VOK1111e9	VOK111	VOK111	1	2	
19	export_B3	2 VOK1111e9	VOK111	VOK111	1	3	
20	export_B4	2 VOK1111e9	VOK111	VOK111	1	4	
21	export_B5	2 VOK1111e9	VOK111	VOK111	1	5	
22	export_B6	2 VOK1111e9	VOK111	VOK111	1	6	
23	export_B7	2 VOK1111e9	VOK111	VOK111	1	7	
24	export_B8	2 VOK1111e9	VOK111	VOK111	1	8	
25	export_B9	2 VOK1111e9	VOK111	VOK111	1	9	
26	export_C10	3 VOK1301e9	VOK130	VOK130	1	10	
27	export_C11	3 VOK1301e9	VOK130	VOK130	1	11	
28	export_C12	3 VOK1301e9	VOK130	VOK130	1	12	
29	export_C1	3 VOK1301e9	VOK130	VOK130	1	1	
30	export_C2	3 VOK1301e9	VOK130	VOK130	1	2	
31	export_C3	3 VOK1301e9	VOK130	VOK130	1	3	
32	export_C4	3 VOK1301e9	VOK130	VOK130	1	4	
33	export_C5	3 VOK1301e9	VOK130	VOK130	1	5	
34	export_C6	3 VOK1301e9	VOK130	VOK130	1	6	
35	export_C7	3 VOK1301e9	VOK130	VOK130	1	7	
36	export_C8	3 VOK1301e9	VOK130	VOK130	1	8	
37	export_C9	3 VOK1301e9	VOK130	VOK130	1	9	

7 Controls Worksheet Headings:

- The first column labeled 'Control_Filename_Prefix' gives the filename of the exported control .csv file. Find the desired control .csv file and write the filename prefix in this column. Note that the control .csv file will also be somewhere in the 'Sample_Filename_Prefix' column in the Sample worksheet, so copying and pasting that filename may be easier.
- The second column labeled 'Sample_Control_ID' is related to the same column in the Sample worksheet. In the Controls worksheet, this column assigns a number to the control. In the Sample worksheet, this corresponding number in the 'Sample_Control_ID' column assigns the control .csv file to the sample .csv file. In this example, there is only one control .csv file which is why you see only one unique 'Sample_Control_ID'.
- The third column labeled 'Control_Name' assigns a name to the control .csv file.

- The fourth column labeled 'Sample_Label_Mix_No' refers to the same information as the same column in the Sample worksheet. **There should be a control for every Label Mix No** which is why there are 14 rows that refer to the same control .csv file in this example.
- The other columns can be filled out as desired, but are not critical for the data analysis.

Control_Filename_Prefix	Sample_Control_ID	Control_Name	Sample_Label_Mix_No	Control_Incubation_Time_With_Antibody	Antibody_Wash_Method	Flow_Cytometer
1		1 Blank Bead	1			
2		1 Blank Bead	2			
3		1 Blank Bead	3			
4		1 Blank Bead	4			
5		1 Blank Bead	5			
6		1 Blank Bead	6			
7		1 Blank Bead	7			
8		1 Blank Bead	8			
9		1 Blank Bead	9			
10		1 Blank Bead	10			
11		1 Blank Bead	11			
12		1 Blank Bead	12			
13		1 Blank Bead	13			
14		1 Blank Bead	14			
15						
16						
17						
18						
19						
20						
21						
22						
23						
24						
25						
26						
27						
28						
29						
30						
31						
32						
33						
34						
35						
36						
37						

8 Labelling Worksheet Headings:

- The first column labeled 'Mix_Number' assigns an ID to the different mixes used in the experiment. This ID is used to determine the label mix in both the Sample worksheet and Controls worksheet under the 'Sample_Label_Mix_No' column.
- The second column labeled 'Import_Column_Number' refers to column of the gated and exported .csv data files that will be imported and used for data analysis. **Ideally, only one channel will be exported from FlowJo and this column can then always be filled with ones.**
- The third column labeled 'Label_Target' assigns a target protein to each label mix.

Mix_Number	Import_Column_Number	Label_Target	Label_Incubated_Concentration	Label_Fluorophore	Label_Isotope	Label_Manufacturer	Label_Catalogue_Number
1	1	1 CD10					
2	2	1 CD13					
3	3	1 CD276					
4	4	1 CD16					
5	5	1 B2M					
6	6	1 CD117					
7	7	1 CD92					
8	8	1 CD103					
9	9	1 CD105					
10	10	1 CD151					
11	11	1 CD309					
12	12	1 CD9					
13	13	1 CD63					
14	14	1 CD81					

9 Beads Worksheet and General Worksheet Headings:

- The Beads worksheet will have the information from the Bead Info Section as discussed in step #3. No changes are necessary for the database to pass QC.
- The General worksheet lists the CSV directory as chosen in step #2. Do not change this worksheet as otherwise the software will not know where the .csv files are.

Import the Database

10 Save the filled out database .xlsx file and return to the MPSPASS software. Click on the 'Import Database' button and choose the database file.

Database QC

11 Once the database has been imported, click on the 'Database QC' button.

If QC is unsuccessful, than the associated error will pop-up in the command window below. The error message will offer an explanation for what failed. The database .xlsx file can be opened again and edited to fix the error. Once the database file has been saved again, the database can be once more be imported and subjected to QC.

12 If QC was successful, the software will allow for the import of the data. Click on the 'Import Data' button and the dataset will be constructed.