

Sep 22, 2021

Preparation of Suppressor tRNA

In 1 collection

Anne Zemella¹, Theresa Richter^{1,2}, Lena Thoring¹, Stefan Kubick¹¹Cell-free and Cell-based Bioproduction, Fraunhofer Institute for Cell Therapy and Immunology (IZI), Branch Bioanalytics and Bioprocesses (IZI-BB), Potsdam, Germany;²University of Potsdam, Potsdam, Germany

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.bqnmvme6

Springer Nature Books

Satyavati Kharde
Springer Nature

ABSTRACT

This is part 3.2 of the "A Combined Cell-Free Protein Synthesis and Fluorescence-Based Approach to Investigate GPCR Binding Properties" collection of protocols:

<https://www.protocols.io/view/a-combined-cell-free-protein-synthesis-and-fluores-bqnmvme6>

Collection Abstract: Fluorescent labeling of de novo synthesized proteins is in particular a valuable tool for functional and structural studies of membrane proteins. In this context, we present two methods for the site-specific fluorescent labeling of difficult-to-express membrane proteins in combination with cell-free protein synthesis. The cell-free protein synthesis system is based on Chinese Hamster Ovary Cells (CHO) since this system contains endogenous membrane structures derived from the endoplasmic reticulum. These so-called microsomes enable a direct integration of membrane proteins into a biological membrane. In this protocol the first part describes the fluorescent labeling by using a precharged tRNA, loaded with a fluorescent amino acid. The second part describes the preparation of a modified aminoacyl-tRNA-synthetase and a suppressor tRNA that are applied to the CHO cell-free system to enable the incorporation of a non-canonical amino acid. The reactive group of the non-canonical amino acid is further coupled to a fluorescent dye. Both methods utilize the amber stop codon suppression technology. The successful fluorescent labeling of the model G protein-coupled receptor adenosine A2A (Adora2a) is analyzed by in-gel-fluorescence, a reporter protein assay, and confocal laser scanning microscopy (CLSM). Moreover, a ligand-dependent conformational change of the fluorescently labeled Adora2a was analyzed by bioluminescence resonance energy transfer (BRET).

For **Introduction and Notes**, please see: <https://www.protocols.io/view/a-combined-cell-free-protein-synthesis-and-fluores-bqnmvme6/guidelines>

DOI

dx.doi.org/10.17504/protocols.io.bqnmvme6

EXTERNAL LINK

https://link.springer.com/protocol/10.1007/978-1-4939-9121-1_4

PROTOCOL CITATION

Anne Zemella, Theresa Richter, Lena Thoring, Stefan Kubick 2021. Preparation of Suppressor tRNA.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.bqnmvme6>

COLLECTIONS ⓘ




A Combined Cell-Free Protein Synthesis and Fluorescence-Based Approach to Investigate GPCR Binding Properties

KEYWORDS

Cell-free protein synthesis, G protein-coupled receptor, Protein modification, Non-canonical amino acids, Amber suppression, Confocal laser scanning microscopy

LICENSE

 This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited









CREATED

Dec 12, 2020

LAST MODIFIED

Sep 22, 2021

OWNERSHIP HISTORY

Dec 12, 2020	 Lenny Teytelman	protocols.io
Jul 05, 2021	 Emma Ganley	protocols.io
Jul 06, 2021	 Emma Ganley	protocols.io
Aug 24, 2021	 Satya K	
Aug 26, 2021	 Satyavati Kharde	Springer Nature
Aug 26, 2021	 Satyavati Kharde	Springer Nature
Sep 22, 2021	 Satya K	
Sep 22, 2021	 Satyavati Kharde	Springer Nature

PROTOCOL INTEGER ID

45493

PARENT PROTOCOLS

Part of collection

[A Combined Cell-Free Protein Synthesis and Fluorescence-Based Approach to Investigate GPCR Binding Properties](#)

2.2 Materials for Preparation of Suppressor tRNA

2.2.1 Generation of PCR Product

1. Vector containing the nucleotide sequence of tRNA^{Tyr}CUA (SupF Gene).
2. tRNA^{Tyr}CUA-specific forward primer (5' CgA gCT CgC CCA CCg gAA TTC 3') and 2'-OMe reverse primer (5' Tgg Tgg Tgg ggg AAg gAT TCg 3').
3. 0.2 ml PCR tubes.
4. PCR cyclor.
5. Taq DNA polymerase.
6. Taq buffer.
7. dNTPs.
8. 25 mM MgCl₂.
9. Agarose gel electrophoresis chambers.
10. Agarose.
11. Rotiphorese 10× TBE buffer.
12. DNA stain.
13. DNA ladder.
14. PCR Purification Kit.

2.2.2 Generation of RNA Transcript

1. T7 RNA Polymerase (f.c. 1 U/μl, Agilent).
2. 5× NTP mix containing 18.75 mM ATP, 18.75 mM CTP, 18.75 mM UTP and 7.5 mM GTP.
3. 5× transcription buffer: 400 mM HEPES-KOH, 0.5 mM Spermidine, 50 mM DTE and 75 mM MgCl₂.
4. DNaseI (1 U per μg plasmid DNA).
5. 10× MOPS buffer: 200 mM MOPS, 50 mM NaOAc, 10 mM EDTA, pH 8.0.
6. MOPS sample buffer: 8% (v/v) formaldehyde, 12 ml formamide, 2.4 ml 10× MOPS buffer, 0.05% (v/v) bromophenol blue to a total volume of 24 ml.

2.2.3 RNA Isolation and Folding

1. TRIzol reagent.
2. High Performance Liquid Chromatography (HPLC) grade Chloroform.
3. HPLC grade Isopropyl.
4. 75% Ethanol.
5. Cooled centrifuge.
6. Nanodrop 2000c.

SAFETY WARNINGS

For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

3.2.1 Generation of PCR Product

- 1 For specific and homogenous 3'-ends of the suppressor tRNA, an additional PCR step before transcription reaction is included. Therefore, the reverse primer contains a 2'-OMe-group to prevent unspecific nucleotides at the 3'-end of the tRNA that can be added by the T7 polymerase during transcription reaction. Amplify the template by pipetting in a PCR tube final concentrations of **[M]1 x Taq Buffer** , **[M]0.2 Milimolar (mM) dNTP mix** , **[M]0.5 Micromolar (μM) forward primer** , **[M]0.5 Micromolar (μM) reverse primer** , **[M]2.5 Milimolar (mM) MgCl₂** , **[M]0.01 ng/μl plasmid** and **[M]0.04 U/μl Taq DNA polymerase** . Fill the reaction with water to a final volume of **250 μl** (see **Note 5**).

3.2.2 Generation of RNA Transcript


6h

2 

Thaw the components for in vitro transcription **On ice** and pipette the reaction at **Room temperature** . Mix

[M]1 x transcription buffer , [M]1 x NTP mix , [M]1 U/μl T7 RNA Polymerase and
[M]8 ng/μl template DNA

3 


Fill the reaction with water to the final volume of  500 μl .

4 

6h

Incubate the reaction for  03:00:00 –  06:00:00 at  500 rpm, 37°C .

5 

Centrifuge the RNA at  12000 x g, 00:01:00 und use the supernatant for the DNaseI treatment (see **Note 6**).

6 

Add  1 U DNase I per 1 μg DNA.

7 

10m

Incubate for  500 rpm, 37°C, 00:10:00 .

3.2.3 RNA Isolation and Folding

1h 41m 15s

8 



Handle the TRIzol and chloroform reagent with care and use a fume hood.

Add a threefold volume of TRIzol to the transcription reaction and mix carefully.

9 

5m

Incubate for  00:05:00 at  Room temperature .

10 

15s

Add  200 μl chloroform for  1 mL TRIzol and mix carefully for  00:00:15 by inverting.

11 

3m

Incubate for  00:03:00 at  Room temperature .

12 

Centrifuge at **12000 x g, 4 °C, 00:15:00** . Isolate the aqueous phase (*see* **Note 7**).

13 

Add **500 µl isopropyl** for **1 mL TRIzol** and mix carefully.

14 

3m

Incubate **Overnight** at **4 °C** .

15 

1h

Centrifuge at **15000 x g** at least for **01:00:00** at **4 °C** and discard the supernatant.

16 

30m

Overlay the pellet with **1 mL 75% ethanol** for **1 mL TRIzol** and incubate for **00:30:00** at **-20 °C** .

17 

Centrifuge at **7500 x g, 4 °C, 00:10:00** . Discard the supernatant and air dry the pellet.

18 

Solve the pellet in water. Measure concentration using a NanoDrop and adjust the concentration to **100 Micromolar (µM)** .

19 

Fold the tRNA by slowly decreasing the temperature from **80 °C** to **25 °C** in a PCR cyclor. The tRNA can be stored at **-80 °C** after shock freezing in liquid nitrogen.