

Aug 27, 2025

## Carver et al, Aged Brain Spatial Profiling - CosMx

 In 1 collection

DOI

[dx.doi.org/10.17504/protocols.io.dm6gp3yb5vzp/v1](https://dx.doi.org/10.17504/protocols.io.dm6gp3yb5vzp/v1)

Chase Carver<sup>1</sup>

<sup>1</sup>Mayo Clinic



Chase Carver

Mayo Clinic

---

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.dm6gp3yb5vzp/v1](https://dx.doi.org/10.17504/protocols.io.dm6gp3yb5vzp/v1)

**Protocol Citation:** Chase Carver 2025. Carver et al, Aged Brain Spatial Profiling - CosMx. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.dm6gp3yb5vzp/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 04, 2023

**Last Modified:** August 27, 2025

**Protocol Integer ID:** 88813

**Keywords:** aged brain white matter, modifiable features of aged brain white matter, aged brain spatial profiling, nanostring cosmx platform, associated microglia, cosmx, cosmx this protocol, spatial molecular imaging, process for spatial molecular imaging

**Funders Acknowledgements:**

NIH Cellular Senescence Network

Grant ID: UG3 CA275669-01




## Abstract

This protocol provides the process for spatial molecular imaging using the Nanostring CosMx platform used in Carver et al., "Senescent- and disease-associated microglia are modifiable features of aged brain white matter".



## CosMx Slide Preparation

- 1 Using FFPE sections, follow the protocol found in the document

 MAN-10159-02\_CosMx\_SMI\_Manu... 2.7MB

## CosMx Instrument run

- 2 Follow the protocol found in the document

 MAN-10161-02-1\_CosMx\_SMI\_Instr... 1.8MB

for imaging of flow cells, FOV selection, and running spatial molecular imaging on the instrument.

The Cell Segmentation Profile used for mouse brain RNA is Configuration B. Cell overlay files are generated with the completion of the CosMx imaging and are used to segment 2D cell areas to be used in conjunction with transcript localization data.

## Protocol references

MAN-10159-02\_CosMx\_SMI\_Manual\_Slide\_Preparation\_Manual