



Quantification of tube topography on micrographs of fluorescent staining of multi-layers cell culture.

Bertrand-David Segard¹

¹[Metcela Inc.]

1 Works for me dx.doi.org/10.17504/protocols.io.bhtaj6ie



Bertrand-David Segard

ABSTRACT

The induction or inhibition of angiogenesis can be assessed in vitro by the differentiation of endothelial cells. Indeed, these cells are capable of forming structures resembling a capillary network in culture. This approach is known as Endothelial Tube Formation Assay (ETFA).

This type of analysis is traditionally performed on a low confluence monolayer of cells to assay the effect of soluble compounds on the formation of tubes. However, the modulation of tube formation by cell-cell interactions requires the co-culture of cells that leads to the creation of a multi-layer micro-tissue. Additionally, the identification of cells in the micro-tissue requires fluorescent staining resulting in discontinuities in the signal.

The present protocol allows to analyze the topography of tubes on micrographs of fluorescent staining of multi-layers cell culture.

The limit of this method resides in the fact that the 3D organization of the micro-tissue is not reconstructed and, thus, a part of the tube network is not analyzed.

DOI

dx.doi.org/10.17504/protocols.io.bhtaj6ie

PROTOCOL CITATION

Bertrand-David Segard 2020. Quantification of tube topography on micrographs of fluorescent staining of multi-layers cell culture.. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bhtaj6ie>

KEYWORDS

co-culture, fluorescent staining, tube topography, quantification

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 23, 2020

LAST MODIFIED

Aug 06, 2020

PROTOCOL INTEGER ID

38466

Software

Install Icy (requires JDK 8+).

1



Icy

Windows/MacOS/Linux

[source](#) by BioImage Analysis Lab (Institut Pasteur)

- 1.1 Download the Icy script "ImageJ background subtraction." (file name: protocolfile-imagej-background-subtraction1)
See: <http://icy.bioimageanalysis.org/protocol/imagej-background-subtraction/>
- 1.2 Install the Icy plugin "Membrane Filter."
See: <http://icy.bioimageanalysis.org/plugin/membrane-filter/>
- 1.3 Install the ImageJ plugin "Luts Macros and Tools Updater."
See: <http://image.bio.methods.free.fr/ImageJ/?Luts-Macros-and-Tools-Updater>
- 1.4 Install the ImageJ plugin "Angiogenesis Analyzer for ImageJ."
See: <http://image.bio.methods.free.fr/ImageJ/?Angiogenesis-Analyzer-for-ImageJ>

Preprocessing

- 2 In Icy, open the image to be analyzed.
 - 2.1 Extract the channel of interest (Channel -> Extract -> Channel x).
Save the new image as TIFF.
 - 2.2 Subtract the background (Tools -> Protocols -> Load: protocolfile-imagej-background-subtraction1).
Settings: Rolling=20.
Save the new image as TIFF.
 - 2.3 Reduce the discontinuities and noise on membranes (Plugin -> Membrane Filter).
Settings: default parameters.
Save the new image as TIFF.

Definition of ROI

- 3 In ImageJ (in Icy), open the image obtained after preprocessing.
 - 3.1 Enhance the contrast (Process -> Enhance Contrast).
Settings: Saturated pixels: 0.0%; Normalize; Equalize histogram.
 - 3.2 Subtract the background (Process -> Subtract Background).
Settings: Rolling ball radius: 5.0 pixels; Sliding paraboloid.
Save the new image as TIFF.

Analysis of tube topography

- 4 In ImageJ (in Icy), open the image obtained after the definition of ROI.
 - 4.1 Start Angiogenesis Analyzer.
Settings: select all options except "Suppress isolated elements."
Threshold values: Minimum object size: 10 pixels; Minimum branch size: 25 pixels; Artifactual loop size: 850 pixels; Isolated element size threshold: 50 pixels; Master segment size threshold: 30 pixels; Iteration number (advised 2 to 5): 3 iterations; Show iteration (for single analysis): 3 iterations.
 - 4.2 Apply the Blurred Mask Tool (Angiogenesis Analyzer -> Blurred Mask Tool) to the whole surface of the image.
 - 4.3 Change the color encoding of the image to RGB color (Image -> Type -> RGB color).
 - 4.4 Run the analysis (Angiogenesis Analyzer -> Analyze HUVEC Phase contrast).

Save and analyze

- 5 Save the result files (Angiogenesis Analyzer -> Save Current Analysis). Measurements are saved as a tab-separated values file (with "xls" extension) and can be further analyzed in a spreadsheet.