

Quick Ligation Protocol (M2200) V.3

New England Biolabs¹

¹New England Biolabs

Version 3

Jun 04, 2021

1

Works for me



Share

dx.doi.org/10.17504/protocols.io.bb2qiqdw

New England Biolabs (NEB)

Tech. support phone: +1(800)632-7799 email: info@neb.com

New England Biolabs Tech Support
New England Biolabs

ABSTRACT

This protocols is to be performed with the Quick Ligation Reaction Buffer. Please see the NEB website for more information.

EXTERNAL LINK

<https://www.neb.com/protocols/0001/01/01/quick-ligation-protocol>

DOI

dx.doi.org/10.17504/protocols.io.bb2qiqdw

EXTERNAL LINK

<https://www.neb.com/protocols/0001/01/01/quick-ligation-protocol>

PROTOCOL CITATION

New England Biolabs 2021. Quick Ligation Protocol (M2200). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bb2qiqdw>

Version created by New England Biolabs Tech Support

KEYWORDS

ligation

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jan 31, 2020

LAST MODIFIED

Jun 04, 2021

OWNERSHIP HISTORY

Jan 31, 2020



New England Biolabs

New England Biolabs

Jun 04, 2021



New England Biolabs Tech Support

New England Biolabs

Jun 04, 2021



New England Biolabs

New England Biolabs

Jun 04, 2021



New England Biolabs Tech Support

New England Biolabs

PROTOCOL INTEGER ID

32560

MATERIALS TEXT

MATERIALS

 [Quick Ligation Kit - 30 rxns](#) **New England**

Biolabs Catalog #M2200S


SAFETY WARNINGS

Please see SDS (Safety Data Sheet) for hazards and safety warnings.

BEFORE STARTING

Thaw and resuspend Quick Ligase Reaction Buffer at  **Room temperature** .

1 

Set up the following reaction in a microcentrifuge tube  **On ice** .

Quick Ligase should be added last. Note that the table shows a ligation using a molar ratio of 1:3 vector to insert for the indicated DNA sizes.

(Use [NEBioCalculator](#) to calculate molar ratios.)

| COMPONENT | 20 µl REACTION |
|------------------------------------|-------------------------|
| Quick Ligase Reaction Buffer (2X)* | 10 µl |
| Vector DNA (3 kb) | 50 ng (0.020 pmol) |
| Insert DNA (1 kb) | 37.5 ng (0.060 pmol) |
| Nuclease-free Water | to 20 µl |
| Quick Ligase | 1 µl |

*The Quick Ligase Reaction Buffer should be thawed and resuspended at room temperature.






2 

Gently mix the reaction by pipetting up and down and microfuge briefly.

3 

Incubate at  **Room temperature** ( **25 °C**) for  **00:05:00** .

4 

Chill  **On ice** and transform  **1 µl** –  **5 µl** of the reaction into  **50 µl** competent cells. Alternatively, Store at  **-20 °C** .

Do not heat inactivate – heat inactivation dramatically reduces transformation efficiency.