



Apr 08, 2022

HuBMAP | GE/ University of Washington Cell DIVE™ Modality Overview

Liz McDonough¹

¹GE Research

1



dx.doi.org/10.17504/protocols.io.3byl4bwmzvo5/v1

Human BioMolecular Atlas Program (HuBMAP) Method Development Community

GE Research



Liz McDonough
GE Research

This is an overview of all protocols currently in use for the GE/University of Washington Cell DIVE collaboration for the Human BioMolecular Atlas Program (HuBMAP). It includes links to each of the individual protocols that make up this project workflow.

DOI

dx.doi.org/10.17504/protocols.io.3byl4bwmzvo5/v1

Liz McDonough 2022. HuBMAP | GE/ University of Washington Cell DIVE™
Modality Overview. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.3byl4bwmzvo5/v1>



Human Biomolecular Atlas Program (HuBMAP)
Grant ID: 3UH3CA246594-02S1

protocol ,

Apr 07, 2022

Apr 08, 2022

60465

1 [Protocol for Tissue Collection from Organ Procurement Organization](#)

- 2 Prepare paraffin blocks and FFPE sections from tissue samples.

[HuBMAP Tissue Preservation Protocol](#)

- 3 Deparaffinize and rehydrate slides.

[Cell DIVE™ Platform | Slide Clearing and Antigen Retrieval](#)

- 4 Characterize antibodies (primary/secondary, direct conjugates, and zenon labelled antibodies) and determine any antigen effects from the Cell DIVE dye inactivation process.

[Cell DIVE™ Platform | Antibody Characterization for Multiplexing](#)

[Cell DIVE™ Platform | Antibody Staining & Imaging](#)

- 5 Prepare direct conjugates for study.

[Cell DIVE™ Platform | Antibody Purification Chemistry](#)

[Cell DIVE™ Platform | Ab Conjugation: Initial Conjugation & Scale up Conjugation](#)

- 6 Perform Cell DIVE™ multiplexed data acquisition on the final cohort.

Staining is done using the Leica Bond MAX and images are acquired on the Leica Cell DIVE imager.