

Jul 10, 2024

Open field test

DOI

dx.doi.org/10.17504/protocols.io.x54v923m4l3e/v1

Chuyu Chen¹

¹Northwestern University, Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815

ASAP Collaborative Rese...

Parisiadou lab



Chuyu Chen

Northwestern University, Aligning Science Across Parkinson's...

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.x54v923m4l3e/v1

Protocol Citation: Chuyu Chen 2024. Open field test. protocols.io <https://dx.doi.org/10.17504/protocols.io.x54v923m4l3e/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 10, 2024

Last Modified: July 10, 2024

Protocol Integer ID: 103170

Keywords: ASAPCRN



Funders Acknowledgement:

**Aligning Science Across
Parkinson's [ASAP-020600]
through the Michael J. Fox
Foundation for Parkinson's
Research (MJFF)**

Grant ID: ASAP-020600

Abstract

The open field test was developed to test locomotor activity.



- 1 Mice were put in a 56 ×56 cm open-field arenas in noise-canceling boxes, illuminated by dim red lights
- 2 The 20-minute-long session started recoding when mice were placed in the center of the arena.
- 3 Locomotor activity was analyzed by the LimeLight 5 (Actimetrics) software and reported as distance traveled.