



JAN 22, 2024

🌐 LB borth or agar

📁 In 1 collection

Adrien Assie¹, Buck Samuel¹

¹Baylor College of Medicine

Samuel Lab



Adrien Assie
Baylor College of Medicine

OPEN  ACCESS



ABSTRACT

LB broth or agar recipe

Protocol Citation: Adrien Assie, Buck Samuel 2024. LB borth or agar. [protocols.io](https://protocols.io/view/lb-borth-or-agar-zy6f7ze)
<https://protocols.io/view/lb-borth-or-agar-zy6f7ze>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Apr 10, 2019

Last Modified: Jan 22, 2024


PROTOCOL integer ID: 22270


1 Start with  900 mL water


2  10 g Bacto-tryptone

3  5 g Bacto-yeast

4  5 g NaCl

5 Adjust pH to 7.0 using  5 Molarity (M) NaOH

6 For LB Agar add  15 g agar

7 Adjust water to  1000 mL

8 Autoclave